

How to Handle Aliquoted Plates of iTru i5 & i7 Primers

Troy Kieran, Todd Pierson, & Travis C. Glenn

contact: Travis C. Glenn

e-mail: travisg@uga.edu

Originally written: November 2013

Last Update: 1 November, 2017

Primer aliquot plates have 1.25 nmoles of dried primer in each well. You will need to reconstitute them to the appropriate volume (125 μ L \rightarrow 10 μ M). The iTru5 sets of 8 are in columns (A-H) and iTru7 sets of 12 are in rows (1-12); refer to the relevant plate map.

Acceptable liquids to use for reconstitution are:

- 1) 10 mM Tris pH 8 (or pH 7.5), or
- 2) TLE (home made or IDTE from IDT = 10 mM Tris pH 8 & 0.1 mM EDTA)

Protocol:

- 1) **Centrifuge the plates to get all the primer to the bottom of the wells.**
- 2) To limit contamination, peel back the foil cover from the plate one row or column (depending on loading scheme) at a time to reconstitute.
- 3) Add 125 μ L of your chosen liquid (e.g., TLE) to each well. Then, using a multichannel pipetter:
 - Skloosh (i.e., pipette up & down) several times & use the pipet tip to help scrape the bottom of the well to dislodge any of the primer that is stuck.
 - Skloosh several more times.
 - Wait a few minutes and then skloosh several more times.
 - Let the primers sit in the liquid at room temperature for at least 5 minutes.
- 4) Once the primers have been fully dissolved, transfer the liquid to new strip tubes.
 - It is wise to use tubes of different color for each set. Be certain to maintain left/right orientation of the strips.
 - It is wise to aliquot each row/column into multiple strips (i.e., so you have <125 μ L per primer in each tube). You want these primers to undergo freezing & thawing as few times as possible, so small aliquots are best (in practice 25 - 62.5 μ L aliquots are reasonable).
- 5) Each primer is now at 10 μ M, and needs to be diluted down to 5 μ M (25 μ L TLE and 25 μ L 10 μ M primer for a 50 μ L aliquot). Now they are ready to use.
 - It is wise to Nanodrop or otherwise quantify your primers. These primer aliquots are quite small, so they are easy to miss. If you miss the tiny dried amount in a well, your primers aren't at 5 μ M & your library PCRs won't work!
 - Use 2.5 μ L for each 50 μ L reaction; you started with enough for 100 reactions.

6) Store the primers at -20°C

7) Take primers out to thaw shortly before use and **skloosh well before using!**