# Quarknet Manual

## CONTENT:

1. Quarknet Experiments
2. Quarknet DAQ卡指令
3. QuarkNet System電腦操作取數據的步驟
4. 解讀QuarkNet數據卡輸出的數據 07/26/10
5. 如何管理QuarkNet e-Lab 帳號
6. TW-QuarkNet-Data-Preparation
7. GPS firmware update procedure
8. 2016 QUARKNET-TW維修工作坊  
維修Quarknet閃爍計數器(Scintillator Counter)  
*拆解, 檢測, 組裝與量測*
9. RP的應用/利用樹莓派的實作IV  
夸克網的應用/硬體篇
10. RP的應用/利用樹莓派的實作V  
夸克網的應用/軟體篇
11. 樹莓派實作(IV, V)線路圖, Linux指令集
12. MACAO-QUARKNET WORKSHOP 2018
13. QuarkNet [Research](#) Workshop 2020, Macau Science Center

## Quarknet Experiments

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TL 4 500  
WC 00 3F

```

[V1
Run Mode      : Off          CE (cnt enable), CD (cnt disable)
Ch(s) Enabled : 3,2,1,0    Cmd DC Reg C0 using (bits 3-0)
Veto Enable   : Off
Veto Select   : Ch0        Cmd DC Reg C0 using (bits 7,6)
Coincidence 1-4: 2-Fold    Cmd DC Reg C0 using (bits 5,4)
Pipe Line Delay: 40 nS     Cmd DT Reg T1=rDelay Reg T2=wDelay 10nS/cnt
Gate Width    : 100 nS     Cmd DC Reg C2=LowByte Reg C3=HighByte 10nS/cnt
Ch0 Threshold : 0.600 vlts
Ch1 Threshold : 0.600 vlts
Ch2 Threshold : 0.600 vlts
Ch3 Threshold : 0.600 vlts
Test Pulser Vlt: 3.000 vlts
Test Pulse Ena : Off

Example line for 1 of 4 channels. (Line Drawing, Not to Scale)
Input Pulse edges (begin/end) set rising/falling tags bits.
-----~----- Input Pulse, Gate cycle begins
-----~----- Delayed Rise Edge 'RE' Tag Bit
-----~----- Delayed Fall Edge 'FE' Tag Bit
-----~----- Tag Bits delayed by PipeLnDly
PipeLineDelay : 40nS
-----|-----|----- Capture Window: 60nS
-----|-----|----- Gate Width : 100nS

If 'RE','FE' are outside Capture Window, data tag bit(s) will be missing.
CaptureWindow = GateWidth - PipeLineDelay
The default Pipe Line Delay is 40nS, default Gate Width is 100nS.
Setup CMD sequence for Pipeline Delay. CD, WT 1 0, WT 2 nn (10nS/cnt)
Setup CMD sequence for Gate Width. CD, WC 2 nn(10nS/cnt), WC 3 nn (2.56uS/cnt)

```

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(National\_Kaohsiung\_Normal\_University\_NKNU)的Lead Teacher

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

# Cosmic Ray e-Lab

Teacher Home Student Home

**High school students use cutting-edge tools to do scientific investigations.**

The Cosmic Ray e-Lab provides an online environment in which students experience the excitement of scientific collaboration in this series of investigations into high-energy cosmic rays. Schools with cosmic ray detectors upload data to a "virtual data grid" portal where ALL the data resides. This approach allows students to analyze a much larger body of


# Cosmic Ray e-Lab

Logged in as group: [CRncku](#) [Logout](#)  
[Helpdesk](#) [My Logbook](#)

Teacher Home Community Standards Site Index Registration  
 Teacher Home Student Home

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# Cosmic Ray e-Lab

Logged in as group: [CRncku](#) [Logout](#)  
[Helpdesk](#) [My Logbook](#)

Teacher Home Community Standards Site Index Registration  
 Teacher Home Student Home

## e-Lab Registration

[Register student research groups with a spreadsheet.](#) Using an existing spreadsheet with your class lists is the easiest way to register groups. The spreadsheet must contain the first and last names of each student in the project and the name of each research group. Decide what role you want the research group to have--user or upload. If your students will be taking the pre- and post-tests for assessment, select **Yes** for the column **In survey**. A template on the next page will get you started.

[Register student research groups.](#) Use this page to register less than ten students. You can link new students to existing research groups or you can create new research groups for students.

[Update your previously created groups.](#) Update your research groups.

[Assign your research groups to new e-Labs.](#) When you register a group, they are registered for the current e-Lab. If you want to register them for new e-Labs as more become available, you need to use this page.

[Update detector IDs for your group.](#) When you analyze data in the Cosmic Ray Elab, most of the analyses need to know information about your detector. This page allows you to assign detector IDs to your group.

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**Update your groups**

Every [DAQ](#) board has a detector id (DAQ number) based on the serial number on the board. Students select it when they upload data. To find it, type "SN" when connected with your detector.

CRncku

Group Name: CRncku  
 Detector ID(s) for CRncku  
 (comma or space separated):

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- Register new students below. To register more students, click on the + button.
- We will create new groups and their associated passwords for you.
- Select  if you want to grant the new group upload permissions for your detectors.
- Select  if you want the new group to take the pretest.

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Your student registration completed successfully.

The groups we created for you (and their associated passwords) are listed below.

If one of the groups you requested already existed in our project, your group name was altered slightly to ensure uniqueness. You may now use the File...Save feature in your browser to save the information below.

Group Name	Password
SSHsiao	knisarro

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Cosmic Ray e-Lab

Logged in as group: [1WQuarkNet](#) [Logout](#) [Helpdesk](#) [My Logbook](#)

Teacher Home Student Home Community Standards Site Index Registration

### Update your groups

SShsiao [Show Group Info](#)


Group Name: SShsiao  
 Academic Year: 2009-2010  
 Role: upload  
 In survey:  No  Yes  
 Password:   
 Verify Password:

Students to delete from "SShsiao":  shsiao1

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Cosmic Ray e-Lab

Logged in as group: [SShsiao](#) [Logout](#) [My Logbook](#)

Home Library Upload Data Posters Site Map Assessment  
 Data Geometry


### Upload raw data collected by your cosmic ray detector.

- Select the detector associated with the data you are uploading.
- Click **Choose File/Browse** to locate the data file on your computer.
- Click **Upload** to upload the file.

Please do not upload files larger than 2 GB in size. You'll have to split them up into smaller pieces. Questions? See the [FAQ](#)

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A **Benchmark File** represents a high quality one day data file for a single geometry and coincidence.

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Other data files will be compared to this file to **bless** the new upload.

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When the detector is reconfigured with a new coincidence or geometry you have to designate a new benchmark (See "Add Benchmark").

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**How do you judge a "Benchmark File"?**

Steady singles & trigger rates.

Sufficient rates to assure adequate S/N.

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Note: for the "6000" series QuarkNet counters, the singles rate is usually between 12 to 35 Hz.

Other style counters will be different and is best determined by plateauing counters (see CR e-Lab→LIBRARY→Resources→Tutorials for plateauing instructions).

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
Data Blessing checks the rates in each channel and the trigger rate for consistency. Blessing involves two stages in your workflow depending on whether you are uploading data with a new configuration or with the same configuration.

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These two stages use two methods for you to "bless" your data and make it available to e-Lab users. Each requires a benchmark file.

- **Manual blessing:** In manual blessing, you have an opportunity to bless files that did not have a benchmark when you uploaded them (Used in stage 1 as shown in the workflow above).
- **Automatic blessing:** In automatic blessing, we compare your upload to a benchmark file you identify (Used in stage 2 as shown in the workflow above).

### Choosing a Benchmark

To select a benchmark file, click on the Benchmark menu item and inspect the blessing plots for several existing uploads. After you inspect the plots, choose an appropriate file to label as benchmark. Read this [tutorial](#) for more information.

You MUST change the benchmark file when you change the configuration (coincidence or geometry) of your detector. You can use the same benchmark

over and over if you never change the coincidence or geometry of your detector.


### Manual Blessing

You can use this method to bless files that you uploaded without automatic blessing. This method still requires you to select a benchmark file or gives you an opportunity to bless data that had been compared to the wrong benchmark at upload time.

Read the tutorial on [manual blessing](#)

### Automatic Blessing

On each data upload you will select a "benchmark" file to compare your data with. You define the benchmark files; they have consistent data rates and similar geometry and triggers to the subsequent upload files that you wish to compare them to.

## Benchmark Tutorial 如何設定腳本

Go to: "Cosmic Ray e-Lab→Upload →Benchmark"

This is where a detector owner will select a **best standard** data file for comparison to other files in order to bless the other. The file owner needs to review the blessing plots and use judgement to select one standard file.

So, first pull down the DAQ list to see your DAQ#s,

<b>Detector</b> Choose detector ▾	<b>Date Range</b> 01/17/2014  to 02/17/2014 	Select Benchmark 
--------------------------------------	--	--


and select your appropriate DAQ#.

Now, set the date range to include one CRMD configuration (coincidence and geometry).

Click "Select Benchmark". This will present a list of candidates for selection.


**Select benchmark files for detector: 6119**

- Navigate through data files to **examine charts** before you choose a benchmark file.
- Select **benchmark** file.
- Add a meaningful **label** to this benchmark file.


Review the list of candidates by clicking on the file number (such as "6119.2013.1012.0") and looking at the Singles Rates and Trigger Rate. Use your judgment to decide if this file is representative of all the data for your particular CRMD configuration. A high quality file will have steady Singles & Trigger Rates, plus Singles will fall between 12 to 35 Hz for the standard QuarkNet counters. Here is an example of good blessing plot:

6202GoodBlessExample

and this shows a bad blessing plot:

6667BadBlessExample

Notice how one of the channels wanders around indicating irregular rates.

After you have selected a representative file, type a helpful name that you will recognize later and click "Save".

Now, when you go back to "Upload" → "Upload" the chosen benchmark will appear on the pull-down list next to the DAQ#.

For the DAQ# for which you are uploading data, select a benchmark file from the pull-down list. This benchmark file will now be compared to the uploaded data to measure its quality.

Later, if the detector configuration changes (coincidence, geometry) then a new benchmark file needs to be selected.

When selecting files for analysis, review the blessing plots by clicking on the "star". See if you agree with the owner assessment for quality data.

## Tutorial for Manual Blessing of Data Files如何將數據檔以手動方式來對角色


To manually "bless" a data file or a set of files, first select the DAQ# from the pull-down list:

**Detector**


**Retrieve all**

To display the list of potential files to be "blessed", click the Retrieve all box. The list will then appear though it may take a moment to generate the list; please be patient.


From the candidates list, you can select individual files or "Bless all files" will attempt to verify the entire list. Shown are two files selected. The singles and trigger rate plots for the highlighted file are displayed.

From the Benchmark pull-down list, select the standard file you wish to compare to the candidate file:

and now click "Check Selected Files" and the CR e-Lab will compare the rates and errors to see if passes.

In this case, one file was "blessed" and the one file failed.

**6420.2013.1029.0** has NOT been blessed. Fail reason: This file failed at: 50550(14:2:30) because the benchmark channel 4 rate: 28.569793282251904 (metadata value) was not between the ranges of comparison set by 14 and 4(10.0 and 18.0) - for these last values, look at the .bless file of the just split file.

**6420.2013.1030.0** has been blessed.

**Hardware notes: 24MHz clock on board, 1.25ns/tick , on board 10nsDT, 50nsDT-GPS**

**Software notes: muonic (Linux and python based), EQUIP(windows,Mac and Java based)**

So, for each of these 2 character hexadecimal words, you essentially have 8 binary bits of information.

The way that this information is encoded is as follows: Bits 0, 1, 2, 3, and 4 denote an even more precise clock count on the board of 1.25 nanoseconds. Bit 5 should be 1 for every valid event, while bit 7 should ONLY be 1 for word 2 - it denotes the start event of a cluster of events. Bits 6 and 7 are always 0. Therefore, the binary word 1011 0010 would mean that a new event in a new cluster has started at 22.5 nanoseconds after the initial time measurement.

Looking at one channel may require you to look between several lines in a cluster. The first line of the

19 cluster is always indicated by the 1 in Word 2, and typically it is easy to follow along the cluster by looking at Word 1 and seeing when the hexadecimal stamps stop being consecutive.

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### Download Firmware Update

1. On a PC, download the update files at <http://www.conwin.com/support.html>
2. The filename is "wi125\_update\_v1.79.00.29.14.rar".
3. Unpack all files to a local working directory. (password: navsync). WinZip may work to unpack, but Dave and Mark downloaded PeaZip.

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### Setup for GPS Card

1. Access GPS card by removing plastic case cover.  
打開GPS外盒
2. Note the location of JP5 (Figure1). Older versions have R9 surface mount pads instead. 參  
考圖一, 找到GPS卡上JP5的所在. 舊版的GPS卡只有R9(表面鍍層)
3. Done with the GPS card for now.  
GPS的初始步驟完成.

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You can remove the short out jumper on the GPS card after "Writing flash memory" appears.

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```

sensor_level = ReadChannel(sensor_channel)
sensor_volts = ConvertVolts(sensor_level,2)

# Print out results
print "-----"
print("Sensor : {} ({}V)".format(sensor_level,sensor_volts))

# Wait before repeating loop
time.sleep(delay)

```

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```

import matplotlib
matplotlib.use("TkAgg")
from pylab import *
import numpy as np
import spidev
import time
import os

# Open SPI bus
spi = spidev.SpiDev()
spi.open(0,0)

# Function to read SPI data from MCP3008 chip
# Channel must be an integer 0-7
def ReadChannel(channel):
    adc = spi.xfer2([1,(8+channel)<<4,0])
    data = ((adc[1]&3) << 8) + adc[2]
    return data

# Function to convert data to voltage level,
# rounded to specified number of decimal places.
def ConvertVolts(data,places):
    volts = (data * 3.3) / 1023
    volts = round(volts,places)
    return volts

# Define sensor channels
light_channel = 0

NT=300
t=np.zeros(NT)
for i in range(NT):

```

```

    t[i]=i
ax=np.zeros(NT)

fig = plt.figure()

bx = fig.add_subplot(111)
plt.subplots_adjust(hspace = .2)
li1, = bx.plot(t, ax,'r-')
plt.grid(True)          #Turn the grid on^M
plt.ylabel('V')
plt.ylim([1.1,1.7])
bx.axes.get_xaxis().set_visible(True)
fig.canvas.draw()
plt.show(block=False)

# Read the light sensor data
cond=True
#while cond:
for k in range(1000):
    tin=time.time()
    for i in range(NT):
        light_level = ReadChannel(light_channel)
        x = ConvertVolts(light_level,2)
        ax[i]=x
# print time.time()-tin
li1.set_data(t,ax)
plt.subplot(111)
fig.canvas.draw()
time.sleep(0.001)

```

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```

import RPi.GPIO as GPIO
import time

GPIO.setmode(GPIO.BCM)

enable_pin = 18
coil_A_1_pin = 4
coil_A_2_pin = 17
coil_B_1_pin = 23
coil_B_2_pin = 24

GPIO.setup(enable_pin, GPIO.OUT)
GPIO.setup(coil_A_1_pin, GPIO.OUT)
GPIO.setup(coil_A_2_pin, GPIO.OUT)

```

```
GPIO.setup(coil_B_1_pin, GPIO.OUT)
GPIO.setup(coil_B_2_pin, GPIO.OUT)
```

```
GPIO.output(enable_pin, 1)
```

```
def forward(delay, steps):
```

```
    for i in range(0, steps):
```

```
        setStep(0, 0, 0, 1)
```

```
        time.sleep(delay)
```

```
        setStep(0, 0, 1, 1)
```

```
        time.sleep(delay)
```

```
        setStep(0, 0, 1, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 1, 1, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 1, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(1, 1, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(1, 0, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(1, 0, 0, 1)
```

```
        time.sleep(delay)
```

```
def backwards(delay, steps):
```

```
    for i in range(0, steps):
```

```
        setStep(1, 0, 0, 1)
```

```
        time.sleep(delay)
```

```
        setStep(1, 0, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(1, 1, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 1, 0, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 1, 1, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 0, 1, 0)
```

```
        time.sleep(delay)
```

```
        setStep(0, 0, 1, 1)
```

```
        time.sleep(delay)
```

```
        setStep(0, 0, 0, 1)
```

```
        time.sleep(delay)
```

```
def setStep(w1, w2, w3, w4):
```

```
    GPIO.output(coil_A_1_pin, w1)
```

```
    GPIO.output(coil_A_2_pin, w2)
```

```
    GPIO.output(coil_B_1_pin, w3)
```

```
    GPIO.output(coil_B_2_pin, w4)
```

```
try:
```

```
    while True:
```

```
        delay = raw_input("Delay between steps (milliseconds)?")
```

```
steps = raw_input("How many steps forward? ")
forward(int(delay) / 1000.0, int(steps))
steps = raw_input("How many steps backwards? ")
backwards(int(delay) / 1000.0, int(steps))
except KeyboardInterrupt:
    setStep(0,0,0,0)
```

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WC 02 0A  
WC 03 00

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H1 Help One  
H2 Help Two  
DG Get GPS Data  
DS Scalar Counts  
DC Get Configuration Information  
DT Time Control  
BA Get Barometer  
TH Get Temperature  
TI Get Time  
V1 View Registers  
V2 View Voltage  
RB Reset Board (counter)  
CE Enable counter to begin data stream (or use Method #1 above) (**Immediately** after RB command)  
(Allow data to stream as long as desired)  
CD Disable counter to stop data stream  
DG Get GPS Data (immediately after CD)  
SA n Save setup, 0=(TMC disable), 1=(TMC enable), 2=(Restore Defaults).  
TL c d Threshold Level, signal ch(0-3)(4=setAll), data(0-4095mV).

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```

import RPi.GPIO as GPIO
import time

```

```

GPIO.setmode(GPIO.BCM)

```

```

1. coil_A_1_pin = 4
2. coil_A_2_pin = 17
3. coil_B_1_pin = 23
4. coil_B_2_pin = 24

```

```

GPIO.setup(coil_A_1_pin, GPIO.OUT)
GPIO.setup(coil_A_2_pin, GPIO.OUT)
GPIO.setup(coil_B_1_pin, GPIO.OUT)
GPIO.setup(coil_B_2_pin, GPIO.OUT)

```

```

def forward(delay, steps):

```

```

    for i in range(0, steps):

```

```

        setStep(0, 0, 0, 1)

```

```

        time.sleep(delay)

```

```

        setStep(0, 0, 1, 1)

```

```

        time.sleep(delay)

```

```

        setStep(0, 0, 1, 0)

```

```

        time.sleep(delay)

```

```

        setStep(0, 1, 1, 0)

```

```

        time.sleep(delay)

```

```

        setStep(0, 1, 0, 0)

```

```

        time.sleep(delay)

```

```

        setStep(1, 1, 0, 0)

```

```

        time.sleep(delay)

```

```

        setStep(1, 0, 0, 0)

```

```

        time.sleep(delay)

```

```

        setStep(1, 0, 0, 1)

```

```

        time.sleep(delay)

```

```

def backwards(delay, steps):

```

```

for i in range(0, steps):
    setStep(1, 0, 0, 1)
    time.sleep(delay)
    setStep(1, 0, 0, 0)
    time.sleep(delay)
    setStep(1, 1, 0, 0)
    time.sleep(delay)
    setStep(0, 1, 0, 0)
    time.sleep(delay)
    setStep(0, 1, 1, 0)
    time.sleep(delay)
    setStep(0, 0, 1, 0)
    time.sleep(delay)
    setStep(0, 0, 1, 1)
    time.sleep(delay)
    setStep(0, 0, 0, 1)
    time.sleep(delay)

```

```

def setStep(w1,w2, w3, w4):
    GPIO.output(coil_A_1_pin, w1)
    GPIO.output(coil_A_2_pin, w2)
    GPIO.output(coil_B_1_pin, w3)
    GPIO.output(coil_B_2_pin, w4)

```

```

#delay = raw_input("Delay between steps (milliseconds)?")

```

```

def turn_voltage(delay,steps):
    if steps>0:
        forward(int(delay) / 1000.0, int(steps))
    else:
        steps=-steps
        backwards(int(delay) / 1000.0, int(steps))

```

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


## QuarkNet Research Workshop 2020, Macau Science Center

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### Schedule:

**15:00-16:00 Cosmic Ray e-Lab, part I:**

communicate DAQ card with Raspberry Pi, measure singles flux rates and double coincidence flux rates for plateau, find the working voltages and threshold levels.

**16:00-16:30 Break**

**16:30-17:30 Cosmic Ray e-Lab, part II:**

save and upload data to e-Lab and use the analyzing tools of cosmic ray e-Lab

**17:30-18:00 Raspberry pi and python programs, part I:**

introduce hardwares including a ADC chip MCP3008, an Op Amp chip LM358 and a stepper motor uln2803.

**18:00-19:00 Raspberry pi and python programs, part II:**

various python programs to perform above experiments in Cosmic Ray e-Lab part I and part II.

**19:00-21:30 Dinner。**