

The [16s rRNA gene](#) is often used to identify different species of bacteria. This is because the sequence is conserved enough that few changes can mean a lot but not conserved completely so there are some changes. This protocol will teach you how to amplify the DNA of a bacterial colony or streak from a plate.

Step 1:

- Setup your [PCR](#) reaction by adding:
(For a 50uL PCR reaction)
 - a. 0.5uL of 10uM 515F primer
 - b. 0.5uL of 10uM 1492R primer
 - c. 25uL 2x Master Mix or 10uL of 2x Master Mix
 - d. Distilled water to 50uL (23uL for 2x Master Mix, 38uL for 5x Master Mix)
 - e. Using an inoculation loop or pipette tip scrape a single colony or a tiny bit of bacteria and mix it into the tube. The amount of bacteria should be visible to the naked eye.
- Use the PCR protocol below

1 x 95C 10 Minutes

30 x

95C 15s

55C 1 minute

68C 2 minute

1x 68C 10 minutes

Step 2:

- Run [Gel Electrophoresis](#) on samples to determine success.
 - a. If you have a gel electrophoresis unit run a gel and see if there is a band of DNA. This is so you know your PCR was successful.
 - i. NOTE: If you do not see a band, try using more of your PCR reaction in the well. Normally 5uL is the standard amount. If you still do not see a band, the PCR amplification was not a success. Retry doubling the amount of bacteria and primers you used.
 - b. Use your PCR purification kit to purify the rDNA that you amplified.
 - c. Send your samples off to a company to sequence along with primers. I generally use the R primer for sequencing. You can submit multiple primers for the same sample just keep them in separate tubes and label all your samples! Usually, one will submit a (F)orward and a (R)everse primer.
 - d. Genewiz is a great company to send your sequence to they will even purify the DNA so you can skip the above step **b**. <https://www.genewiz.com/>
- Once they receive the samples they will send you back two files for each sequencing run.

- a. One is usually labeled <name of sample>.seq and <name of sample>.ab1
- b. The .seq file is the one that contains the DNA sequence that they acquired from the run.
- c. The .ab1 file contains a spectrogram.
- d. To look at the ab1 file you need special software, FinchTV is probably one of the more popular ones (<http://www.geospiza.com/Products/finchtv.shtml>).

Step 3:

- In order to identify the bacteria that our sequence belongs to and so identify our bacteria of interest we need to compare our sequence(s) to the database at NCBI (http://blast.ncbi.nlm.nih.gov/Blast.cgi?PROGRAM=blastn&PAGE_TYPE=BlastSearch&LINK_LOC=blasthome)
- Where it says "Enter accession number(s), gi(s) or FASTA sequence(s)" paste your 16s rDNA sequence into that box.
 - a. If the sequence in the file is really short (less than 30 bases) usually it means the sequencing failed.
- Where it says "organism" type in "bacteria". Then go to the bottom and click the button that says "BLAST".
- After 1-30 seconds your results should come up. It is often that the #1 hit is the genus of bacteria. The species can be a little more complicated as generally you won't find a 100% match. If you also did a forward primer sequencing reaction BLAST that one also and compare the species. If both the forward and reverse primers match species there is a good chance that your bacteria belongs to that species or a very closely related one that has not been sequenced and identified yet!