

Photochemical Experiment Product Yield Study: α -Pinene + OH

It is recommended that the following procedure be carried out by each participating chamber as part of the ICARUS inter-comparison study. The purpose of this study is to compare SOA and individual product yields using vapor/particle wall loss corrections and the photon flux protocol to determine modeled vs. experimental OH (via decay of VOC).

I. Suggested Procedures

- α -Pinene photooxidation experiments will be carried out at ambient temperature (~ 295 K), atmospheric pressure (~ 1 atm), and low relative humidity (15% RH). At Caltech, temperature and RH will be monitored with a Vaisala HMM211 probe.
- Experiments will be performed in the presence of $(\text{NH}_4)_2\text{SO}_4$ seed aerosol, and at an initial NO_x mixing ratio of 15 ppb (injected as pure NO). At Caltech, O_3 and NO_x mixing ratios are quantified by a Horiba APOA-360 O_3 absorption monitor and a Teledyne T200 NO_x monitor, respectively. The detection limits for O_3 , NO, and NO_2 are 0.5, 0.4, and 0.4 ppb, respectively.
- Prior to each experiment, the chamber should be flushed with dry, purified air for 24 h such that the particle number and volume concentrations are $< 10 \text{ cm}^{-3}$ and $0.01 \mu\text{m}^3 \text{ cm}^{-3}$, respectively.
- α -Pinene ($\geq 99\%$, Sigma-Aldrich) should be injected into the chamber such that the initial mixing ratio is ~ 100 ppb.
- Polydisperse seed aerosol should be generated via atomization of a dilute (0.06 M) aqueous solution of $(\text{NH}_4)_2\text{SO}_4$, followed by diffusive drying and neutralization. The target volume concentration and mean particle diameter are $\sim 70 \mu\text{m}^3 \text{ cm}^{-3}$ and ~ 100 nm, respectively.
- H_2O_2 should serve as the source of OH on photolysis under broadband UV irradiation. At the slow H_2O_2 photolysis rate in the Caltech chamber (350 nm backlights) and the relatively high initial mixing ratios of H_2O_2 , steady-state OH and HO₂ concentrations are achieved (Eddingsaas 2012). The initial H_2O_2 mixing ratio should produce steady-state OH and HO₂ concentrations of $\sim 2 \times 10^6$ and $\sim 1 \times 10^{10}$ molecules cm^{-3} , respectively. At Caltech, this requires an initial H_2O_2 mixing ratio of ~ 2 ppm.
- Experiments should be run for 4-h.

II. Instrumentation

1. Required Instruments

- An instrument to monitor α -pinene mixing ratio, such as PTR-MS, SIFT-MS, GC-FID, GC-ECD, IMS, or similar.
- An instrument to monitor gaseous oxidation products, such as negative ion CIMS (any reagent ion), PTR-MS, SIFT-MS, or similar.

- An instrument to monitor particle size and concentration, such as SMPS, AMS, or similar.

2. Optional Instruments

- An instrument to monitor particle composition, PILS coupled with HPLC, FIGAERO, or similar.

III. Example Calibration and Data Analysis Protocols for Caltech Chamber

1. Gas Phase α -pinene Mixing Ratio using GC-FID

α -Pinene mixing ratios are quantified with an Agilent 6890N gas chromatograph equipped with a flame ionization detector (GC-FID) and operated with an Agilent HP-5 column (30 m $\times$ 0.32 mm, 0.25 μ m). The GC-FID is calibrated with a commercial α -pinene standard (>99%, Sigma- Aldrich) over a mixing ratio range from 100 to 200 ppb using a gas-tight volumetric syringe and a mass-controlled dilution flow of N₂ into a 100 L Teflon bag.

At Caltech, the GC-FID is also calibrated with ppm-level bags (10-20 ppm) prepared via an analogous method and cross-calibrated using Fourier transform infrared spectroscopy (FT-IR) with tabulated absorption cross sections for α -pinene (Sharpe 2002).

2. SMPS Data Analysis

At Caltech, SMPS data analysis is carried out through data inversion based on finite element analysis of particle trajectories along the TSI 3081A long column DMA, Monte Carlo simulation of particle Brownian motion, and residence time distribution modeling of the TSI 3010 CPC (Mai 2018). Based on instrumental limitations and inversion accuracies, the SMPS used in the Caltech chamber can measure particle size distributions from 20 nm up to 800 nm, which is sufficient to cover the SOA size range in this experiment. Noted that using different DMA transfer functions (Stolzenburg 1988, Russell 1995, Collins 2004, Mamakos 2008, Dubey 2011) and instrumentation setup may result in differences in SOA yield calculations.

As a result of particle growth driven by gas/particle-phase chemistry and gas-particle partitioning, mean geometric diameters typically shift from $\sim$ 100 nm to $\sim$ 230 nm after a 4-hr experiment. Aerosol volume concentrations are calculated from particle size distributions assuming homogenous spherical particles. Aerosol mass concentrations are derived assuming an effective density for α -pinene SOA of 1.25 g mL⁻¹.

3. AMS Data Analysis

AMS ionization efficiency is calibrated using dry, 350 nm NH₄NO₃ particles, generated from a dilute (0.01 M) aqueous solution of NH₄NO₃ and size-selected with a DMA. Detection limits for each class of chemical constituents are calculated as three times the

standard deviation of blank signals measured from HEPA filter samples (~30 min) taken before each experiment.

AMS data are analyzed using the SQUIRREL and PIKA modules for Igor Pro, and are corrected for gas-phase interferences (Allan 2005, Aiken 2008) and composition-dependent collection efficiencies (Middlebrook 2012). Elemental O:C and H:C ratios are calculated using the “Improved-Ambient” elemental analysis method for AMS spectra (Canagaratna 2015).

IV. References:

- Aiken, et al., *O/C and OM/OC Ratios of Primary, Secondary, and Ambient Organic Aerosols with High-Resolution Time-of-Flight Aerosol Mass Spectrometry*, Environmental Science and Technology, 2008, **42**, p. 4478-4485.
- Allan, et al., *A generalised method for the extraction of chemically resolved mass spectra from Aerodyne aerosol mass spectrometer data*, Journal of Aerosol Science, 2005, **35**, p. 909-922.
- Canagaratna et al., *Elemental ratio measurements of organic compounds using aerosol mass spectrometry: characterization, improved calibration, and implications*, Atmospheric Chemistry and Physics, 2015, **15**, p. 253-272.
- Collins et al., *The scanning DMA transfer function*, Aerosol Science and Technology, 2004, **38** (8), p. 833-850.
- Dubey et al., *A new approach to calculate diffusional transfer functions of scanning DMAs*, Aerosol Science and Technology, 2011, **45** (8), p. 1031-1040.
- Eddingsaas et al., *α -pinene photooxidation under controlled chemical conditions – Part 1: Gas phase composition in low- and high- NO_x environments*, 2012, **12**, p. 6489-6504.
- Kenseth et al., *Synergistic $\text{O}_3 + \text{OH}$ oxidation pathway to extremely low-volatility dimers revealed in β -pinene secondary organic aerosol*, Proceedings of the National Academy of Sciences, 2018, **115** (33) p. 8301-8306
- Mai et al., *Scanning DMA Data Analysis I. Classification Transfer Function*, Aerosol Science and Technology, 2018.
- Mai et al., *Scanning DMA data analysis II. Integrated DMA-CPC instrument response and data inversion*, Aerosol Science and Technology, 2018
- Mamakos et al., *Differential mobility analyser transfer functions in scanning mode*, Journal of Aerosol Science, 2008, 39 (3), p. 227-243.
- Middlebrook et al., *Evaluation of Composition-Dependent Collection Efficiencies for the Aerodyne Aerosol Mass Spectrometer using Field Data*, Aerosol Science and Technology, 2012, **46**, p. 258-271.
- Russell et al., *Asymmetric instrument response resulting from mixing effects in accelerated DMA-CPC measurements*, Aerosol Science and Technology, 1995, **23** (4), p. 491-509.

Sharpe et al., *The PNNL quantitative IR database for infrared remote sensing and hyperspectral imaging*, Proceedings of the 31st Applied Imagery Pattern Recognition Workshop, 2002.
<https://secure2.pnl.gov/nsd/nsd.nsf/Welcome>

Stolzenburg et al., *An ultrafine aerosol size distribution measuring system*, Ph.D. Thesis, University of Minnesota, 1988.