



Spence Laboratory - University of Michigan
7/2/2018

Human Tissue Digestion (for downstream FACS isolation of ECs)

Digestion of tissue to single cell state

1. Mince tissue to $<1 \text{ mm}^3$. Incubate in 50 mL conical with **9 mL collagenase II** and **1 mL dispase** (per gram of tissue). See below for reagent prep.
 - Incubate at 37°C for 30 minutes - **continuous** agitation
2. Add **75 μL DNase** (per 10 mL of digestion solution)
 - Incubate at 37°C for 30 minutes - **continuous** agitation
3. Add **5 mL EC isolation medium** (per 10 mL of digestion solution). Pass cells through a **100 μm cell strainer** to remove undigested tissue. Optional rinse with an additional 5 mL EC isolation medium. Collect in a clear 50 mL conical.
 - Centrifuge at 400g for 5 minutes (RT)

Removal of RBCs and dead cells

4. Resuspend cells in **10 mL Ficoll cell separation medium** (per gram starting material).
5. Layer suspension on **7.5 mL Ficoll-Paque PLUS** (prewarmed to RT). To layer, carefully tilt conical on its side and *slowly* pipette the cell suspension down the side of the tube. Centrifuge at 400g for 20 minutes at 18°C **with as slow of acceleration as possible and without brakes**...this will end up taking approximately 45 minutes in total.
 - Do not disturb the Ficoll-Paque - will ruin density gradient
 - Put cells on **ice** immediately after centrifugation
6. Transfer interphase to a fresh tube and add **2 volumes of EC isolation medium**. Remove top layer to more easily access the middle layer containing cells of interest.
 - Centrifuge at 400g for 5 minutes
7. Resuspend the cells in an appropriate medium depending on downstream application.

EC Medium	EC Isolation Medium
• 78 mL RPMI-1640 • 20 mL FBS • 1 mL penicillin-streptomycin (10,000 U/mL; 10,000 mg/mL) • 1 mL L-Gln (200mM - not included if RPMI contains glutamine already) *Filter through a 0.2 µM filter *Store at 4°C	<u>Component A</u> (80% vol/vol) • 74 mL RPMI-1640 • 1 mL L-Gln (200 mM - not included if RPMI contains glutamine already) • 1 mL penicillin-streptomycin (10,000 U/mL; 10,000 mg/mL) • 1 mL sodium pyruvate (100 mM) • 1 mL nonessential amino acids (10 mM) • 1 mL HEPES (1 M) • 1 mL MEM vitamins (100x) • 20 mL FBS
<u>Ficoll-Paque Separation Medium</u> • <u>EC Medium</u> ○ 1% (vol/vol) 2 M sodium citrate	<u>Component B</u> (20% vol/vol) • <u>EC Medium</u> *Store at 4°C for 1 week

- **Collagenase Type II Solution:** 0.1% (w/v) in PBS, Filter sterilized. Made fresh each time.
- **Dispase Solution:** Dissolve dispase to 2.5 units/mL in PBS. Filter sterilized. Made fresh each time.

Protocol adapted from van Beijnum, J. R., Rousch, M., Castermans, K., van der Linden, E., & Griffioen, A. W. (2008). Isolation of endothelial cells from fresh tissues. *Nature Protocols*, 3(6), 1085–1091.
<http://doi.org/10.1038/nprot.2008.71>