

Test Plan
Energy Efficiency
Phase 4, Medical Concentrator Assessment
Draft A

The below tests are to be performed on O2 concentrator models supplied by DT Global specifically for Performance Testing.

Test Plan
Energy Efficiency
Phase 4, Medical Concentrator Assessment
Draft A

1. Test Name

- **Usability Assessment**

2. Relationship with Operation in Low Resource Settings (LRS) – for Reference Only

- Oxygen concentrators are often used around the clock, and can have a large power demand.
- In LRS, units may be operated on solar power.

3. Test Objective

- To assess energy efficiency characteristics of medical oxygen concentrators.

4. Test Apparatus

- Medical oxygen concentrator to be assessed
- Multimeter

Test Plan
Energy Efficiency
Phase 4, Medical Concentrator Assessment
Draft A

5. Test Procedures

5.1 Variable power

1. Criterion: It is advantageous for the unit's power consumption to scale with delivered flow rate.
2. Test method:
 - a. Turn the unit on.
 - b. Using the multimeter, measure the power consumed by the unit at 0, 2, and 10 L/min.
3. Impact: Major preferred. Units that meet the criterion receive a Major Advantage, which may impact the rating of the unit.

5.2 Soft start

1. Criterion: It is advantageous for the unit to have a soft-start feature to avoid high current at startup.
2. Test method:
 - a. Using the multimeter, record the current draw when the unit is powered on.
 - b. We also spoke to the vendor to inquire if the unit had a soft-start feature.
3. Impact: Major preferred. Units that meet the criterion receive a Major Advantage, which may impact the rating of the unit.

5.3 Power consumption at max flow

1. Criterion: It is advantageous for the unit to use less than 400 W at 10 L/min.
2. Test method:
 - a. Turn the unit on and set the flow to 10 L/min.
 - b. Using the multimeter, measure the power consumed.
3. Impact: Minor preferred. Units that meet the criterion receive a Minor Advantage, which is noted in the report but does not affect the rating of the unit.

5.4 Power consumption at low flow

4. Criterion: It is advantageous for the unit to use less than 300 W at 1-5 L/min.
5. Test method:
 - c. Turn the unit on and set the flow to 1 L/min.
 - d. Using the multimeter, measure the power consumed.
 - e. Repeat for 2, 3, 4, and 5 L/min.
6. Impact: Minor preferred. Units that meet the criterion receive a Minor Advantage, which is noted in the report but does not affect the rating of the unit.