

# ALICE Workflow

## Imaging, File Renaming, Cropping & Scale Bar Stamping

**Author:** Louise Allan

**Contributors:** Robyn Crowther, Elizabeth Devenish, Phaedra Kokkini, Nicola Lowndes, Peter Wing

**Institution:** Natural History Museum, London

**Date created:** 28 June 2019

**Last updated:** 2 March 2020

## Image Capture

### 1. Turn on imaging equipment


- Switch on the top LED light (switch located at the top of the lightbox)
- Plug in the additional bottom LED lights
- Turn on all 6 cameras

### 2. Software




- Open the DSLR Multi Camera software:
- Check all 6 cameras have connected
  - check status in the programme window or number of cameras in drop-down list - example 1
- check the camera settings by selecting each camera - example 1 (camera set to manual)
  - ALICE1** = F20; ISO 400; ss 1/80; **ALICE2** = F14; ISO 400; ss 1/50; **ALICE3-6** = F20; ISO 400; ss 1/50
- Open live view for either a single camera (Live View) or all 6 cameras (Multi Live View) - example 2

Example 1:



Example 2:



### 3. Image the specimens

At the start of each session, the software creates a date folder within the '1\_ALICE image capture' folder in the folder 'processing' (only after the first image is taken). All subsequent images taken in a session are saved in this folder.

**Note:** To take an image, click 'Release' or use the spacebar


- Calibration:** Before any specimen images are taken, image the calibration square (these will be needed for any future label transformation)
- Specimens:**
  - insert the corresponding location and taxon IRN label(s) for the specimen in the top left hand corner of the ALICE1 image (raised foam platform in template)
  - Place the specimen in the centralised hole
    - use the right hole for specimens where the pin is through the side of the label(s)
    - use the left hole for specimens where the pin is through the centre of the label(s)

**Note:** Position folder window so the set of six images align with the next six to identify irregularities.

---

## Daily File Processing

There are four main folders that will be used whilst processing (shortcuts on the right hand side of the desktop screen):



### At the end of the day:

#### 1. Calibration Images

##### 1.1 File transfer

- a. **MOVE (cut and paste)** the 6 calibration images from '1\_ALICE image capture' folder to the ALICE folder on Unified - shortcut on the desktop:

\* Folder structure: **desktop \ DCP Unified ALICE folder \ "project" \ YYYY\_MM\_DD \ calibration**

#### 2. Rename Specimen Image Files via BarcodeFiler

##### ● ALICE1 (dorsal):

##### 2.1. File Transfer

- a. Separate the dorsal images from the lateral and label images:
  - i. Within the day's session folder in '1\_ALICE image capture' use the folder search function to lookup 'ALICE1'
  - ii. **MOVE** the ALICE1 dorsal images from '1\_ALICE image capture' folder to the 'input' folder in '2\_BarcodeFiler' (desktop shortcut):

\* Folder structure: **desktop \ 2\_BarcodeFiler \ input \**

##### 2.2. BarcodeFiler



- a. Using the BarcodeFiler software rename the ALICE1 (dorsal) images
  - click 'Run'(use the .ini settings = 'ALICE\_Entom NHMUK\_NO memory between images\_UID\_LOC\_TAX')
  - leave this to run overnight - *continue the renaming the following day (see 2.4 and 2.5)*

##### ● ALICE2-6 (lateral & label):

##### 2.3. File transfer

- a. **COPY** 'ALICE2'- 'ALICE6' (lateral & label) images from '1\_ALICE image capture' to the 'ALICE2-6' folder in '4\_ALICErenaming' (desktop shortcut):

\* Folder structure: **desktop \ 4\_ALICErenaming \ ALICE2-6 \**

**Note:** ensure that the "ALICE2-6" folder is empty before copying the new set of images

**Note:** copy the files rather than moving them so there is a backup if the renaming needs to be repeated

---

## Following day:

- **ALICE1 (dorsal):**

### 2.4. BardecodeFiler

#### a. 'renamed' folder:

The BardecodeFiler software will transfer the renamed files to the BardecodeFiler 'renamed' folder:

\* Folder structure: `desktop \ 2_BardecodeFiler \ renamed \`

- image files will be renamed with the original filename followed by UID\_LocIRN\_TaxIRN

i.e. 0001\_ALICE1-UID\_LocIRN\_TaxIRN

#### Files with issues:

- image files missing one or more of the three barcode values will be renamed with the original filename followed by the barcodes that could be read
- check the image files for short filenames:
  - most common error is UID not read (obscured by specimen / not found)  
i.e. 0001\_ALICE1-\_LocIRN\_TaxIRN
  - manually add the missing information to the filenames of these images  
(use the hand barcode scanner to prevent manual renaming errors)

#### b. 'Exceptions' folder

Check BardecodeFiler to see if any exception images have been created, these can be found in:

`desktop \ 2_BardecodeFiler \ exceptions \`


- if there are images:
  - rename them manually and move them to the 'renamed' folder

Files Created	Exceptions
00365	00014

**Note:** if we are not using the location and taxon temporary IRN label then the filename must include the **double underscores**, otherwise the renaming script and EMu ingestion script will not work:  
i.e. 0001\_ALICE1-UID\_\_

#### c. Delete both the 'barcodefiler' text file, and 'results' Excel file

**Note:** these may affect the renaming script if copied over with the images



#### d. Once the renaming is complete **COPY** all the dorsal images from 'renamed' to the 'ALICE1-Dorsal' folder in '4\_ALICErenaming' (desktop shortcut):

\* Folder structure: `desktop \ 4_ALICErenaming \ ALICE1-Dorsal \`

**Note:** ensure that the "ALICE1-Dorsal" folder is empty before copying the new set of images

**Note:** copy the files rather than moving them so there is a backup of the initial renamed ALICE1 images.

#### e. Barcodefiler may now be closed

- **ALICE 1 (dorsal) and ALICE2-6 (lateral & label):**

### 2.5. Batch Renaming - Script

#### a. Ensure the images are in the correct folders:

- 'ALICE1-Dorsal' = dorsal images
- 'ALICE2-6' = lateral and label images

#### b. Renaming

- double click the "ALICE-rename\_script" file (.bat) (use latest version)
  - script location: `desktop \ 4_ALICErenaming \`

#### c. All files will be renamed = see script explained (next page)

*Script explained:*

using a renaming script 1) the lateral and label images will be renamed using the corresponding dorsal image as a reference, and 2) the dorsal image will have the original file name removed:

File	Start	End
ALICE1	0001_ALICE1-UID_LocIRN_TaxIRN	UID_LocIRN_TaxIRN
ALICE2	0001_ALICE2	UID_lateral
ALICE3	0001_ALICE3	UID_additional_1
ALICE4	0001_ALICE4	UID_additional_2
ALICE5	0001_ALICE5	UID_additional_3
ALICE6	0001_ALICE6	UID_additional_4

## • Quality Checks

### a. Dorsal filenames

check for the following:

- any suffixes i.e. \_1 or (1)
- misread or duplicate UID barcode values

### b. Lateral and label filenames


check the files for the following and correct any issues:

- files that have not been renamed
  - rename these manually using original images for verification (dorsal in 'renamed' folder)
- unusually suffixes i.e. '\_additional\_5', 'additional\_6' etc

### Unusual pinned specimens: multi-barcode specimens require special renaming

Some multi-specimen pins may be assigned multiple barcodes. The specimens will need to be imaged three times, and as a result will show up in both the **dorsal** and **lateral and label** folders as misread and will need to be manually renamed.

**Note:** pay attention to which plane each 'additional\_#' is assigned to when renaming via the table above (eg: The direction of view of 'additional\_1' in one specimen should equal the same direction of view in another).



## 3. File Processing

### Dorsal images:

- if the location and taxon temporary IRN label was present in the initial dorsal image then these images must be cropped - see Stage 3.1.
- if no temporary labels in the initial dorsal image i.e. only the specimen, then no cropping is required:
  - MOVE** the dorsal images to the 'final' folder:
    - \* Folder structure: desktop \ final \ YYYY\_MM\_DD \ images

### Lateral and label images:


- if you wish to add a scale bar to the lateral images - see Stage 3.2.
- if you do not wish to add a scale bar:
  - MOVE** the lateral and label images to the 'final\_ReverseSide' folder:
    - \* Folder structure: desktop \ final\_ReverseSide \ YYYY\_MM\_DD \ images

### 3.1. Cropping - Dorsal Images

#### 3.1.1. XnConvert



- a. For the dorsal images in the 'ALICE1-Dorsal' folder crop out the location and taxon IRN label from the final image using XnConvert:
  - i. drag and drop the images from the folder into the XnConvert 'input' tab
  - ii. select the appropriate project settings from the dropdown list (bottom left - next to the folder symbol)
  - iii. click 'Convert' to crop the images



- a pop up window will appear asking if you wish to replace the files, select "yes to all"

**Note:** some images may need different requirements for cropping - these images can be manually cropped by dragging the individual files to Input, and choosing the required settings in Actions. Finish these images before batch-cropping the rest of the files.

- b. Once cropped **MOVE** the dorsal images to the 'final' folder:
  - \* Folder structure: **desktop \ final \ YYYY\_MM\_DD \ images**

Dorsal image processing is now **complete**.


### 3.2. Scale Bar - Lateral Images

#### 3.2.1. File transfer

- a. Separate the lateral and label images:
  - i. in the folder 'ALICE2-6' use the folder search function to select 'additional'
  - ii. **select** and **MOVE** all the label images to the 'final\_ReverseSide' folder:
    - i.e. UID\_additional\_1, UID\_additional\_2, UID\_additional\_3, UID\_additional\_4
    - \* Folder structure: **desktop \ final\_ReverseSide \ YYYY\_MM\_DD \ images**
- b. You should be **left with** just the lateral images in the 'ALICE2-6' folder
  - i.e. UID\_lateral



#### 3.2.2. ImageJ

- a. Using ImageJ stamp a 5 mm scale bar on the lateral images using a macro script
  - i. open: Plugins → Macros → Run
    - "Run Macro or Script" - select the macro 'ALICE\_800D 60mm\_scale bar\_5mm\_Lateral'
      - \* Folder structure: **SBILPictures \ documents \ ImageJ \ macros**
    - use the  to find the documents folder (3 clicks)
    - "Choose the input folder!" - select the images folder 'ALICE2-6'
      - choose the desktop option to find the folder
    - "Choose the output folder!" - select where to save the new images ('ALICE2-6' folder)
 

**Note:** by saving to the same folder the macro will overwrite the original images
    - The macro will begin running
 

**Note:** there is no indication when the macro has completed - check the final image (help document [Add scale bar using ImageJ](#))
- b. **MOVE** the lateral images to the 'final\_ReverseSide' folder:
  - \* Folder structure: **desktop \ final\_ReverseSide \ YYYY\_MM\_DD \ images**

Lateral and label image processing is now **complete**.

---

## Daily File Cleardown

- a. Once the image renaming and processing is completed **delete** the following copies:
  - i. '1\_Alice image capture' folder area
    - original lateral and label images (also delete the capture date folder)
  - ii. '2\_BardecodeFiler' folder area
    - 'original\_processed' folder = original (non-renamed) images
    - 'renamed' folder = dorsal renamed ALICE1 images
    - 'exceptions' folder = dorsal renamed ALICE1 images

---

## File Transfer

### 1. FLOWS

- a. **Copy** the date folders with the images daily to the corresponding EMu import shared folder - daily:
  - i. **"final"** copy to `\\dfs-ctdb\emu-import-Slides_SLR_X`
  - ii. **"final\_ReverseSide"** copy to `\\dfs-ctdb\emu-import-Slides_SLR_X_Reverse`

**Note:** the Data Manager will be notified about the images and, once in EMu, the FLOWS script will automatically clear down the images and date folders in **"final"** and **"final\_ReverseSide"**

### 2. Future Label Extraction and Transformation (label images)

- a. **Copy** the lateral and label images (additional\_1 to \_4) to the ALICE folder on Unified where the calibration images were saved - shortcut on the desktop:
  - \* Folder structure: `desktop \ DCP Unified ALICE folder \ "project" \ YYYY_MM_DD \`