Investigate the medicinal use of plants by Australian Aboriginal and Torres Strait Islander Peoples in infection control

Name:	Date:							
Aim/Question What do you want to find out?	To investigate the medicinal use of plants by Australian Aboriginal and Torres Strait Islander Peoples in controlling bacterial infection							
Prediction INDIVIDUAL	I predict the agent that will inhibit bacterial growth the best wil ☐ Tea Tree oil ☐ Eucalyptus oil ☐ Manuka honey ☐ Salt water ☐ Rubbing Alcohol	ll be:						
Independent variable What will you change?								
Dependent variable What will you measure?	I will measure:	I will measure this using:						
Controlled variables What will you keep the same?								
Hypothesis	Ifthen							


Non-pathogenic bacteria Permanent marker Gloves **Materials** petri dishes & nutrient Agar Sticky tape Bunsen burner (list all the cotton buds Incubator Matches Scissors equipment you will need - use the Tea Tree oil variables to help you) Filter paper Salt water Heat proof mat Eucalyptus oil **Tweezers** Manuka honey Rubbing Alcohol Risk **Risk** How will I manage the risk assessment What risks can you identify in this investigation? What measures can be taken to minimise the risks?

Method

1. Collect your equipment and set up the Bunsen Burner.

NOTE: Do not remove the lid of the Petri dish until you are ready to inoculate it

- 2. On the lid of the Petri dish write your name, group members names and the sample you are testing around the edge of the lid
- 3. Using a sterile cotton bud collect some bacterial sample, remove the lid of the Petri dish and LIGHTLY rub the cotton bud on the agar in a zigzag pattern across the whole plate and put the lid back on. Dispose of the cotton bud.
- 4. IF YOU ARE MAKING THE CONTROL GO TO STEP 8
- 5. Light the Bunsen Burner and turn to blue flame.
- 6. Wave the tweezers in the flame to kill any bacteria that may be on them, DO NOT overheat them
- 7. Pick up a piece of filter paper that has been soaking in your sample and place it at 12 o'clock, repeat this so there is a piece of sample at 3 o'clock, 6 o'clock and 9 o'clock. Put the lid back on the Petri dish.



- 8. Wrap sticky tape around the edge of the Petri dish securing the lid to the base to prevent the lid from being removed again.
- Place the Petri dish in the incubator with the lid down to prevent any water droplets from falling on the agar and culture
- 10. Leave for 24 48 hours
- 11. Record observations and the number of colonies growing on your dish. Share the results with the class. Take photos of the results.

Data collection Antibacterial **Number of colonies TABLE** agent Test 1 Test 2 Test 3 Test 4 **Average** Control Tea tree oil **Eucalyptus** Oil Manuka honey Salt water Rubbing Alcohol **Observations** Antibacterial Diagram / Photo **Observations** agent Control Tea tree Oil **Eucalyptus Oil** Manuka honey Salt water Rubbing Alcohol

GRAPH																						
Graph your		Ш																				
results using a line graph.		Н			Н	\dashv				\vdash	\dashv		-	_				H	Н	Н	\dashv	
		Н			Н	\dashv				\dashv	\dashv		\dashv	_					Н	Н	\dashv	
The vertical axis		Н								\forall	\dashv						\vdash		Н	Н	\dashv	
is the average number of																						
colonies.																						
The horizontal axis is the		Н				-				Н	\dashv		\Box	_			L	L	Н	Н	\dashv	
antibacterial		Н				\dashv				\dashv	\dashv		-	_			H	H	Н	Н	\dashv	
agent		Н				\dashv			Н	\dashv	\dashv						\vdash		Н	Н	\dashv	
		П								\Box	\exists									П	\dashv	
INDIVIDUAL		Ш																	Ш	Ш		
		Н				\dashv				\dashv	\dashv		-	_		_	H	H	Н	Н	\dashv	
		Н			Н	\dashv				\dashv	\dashv		\dashv						Н	Н	\dashv	
		Н								\forall	\exists						\vdash		Н	П	\dashv	
Discussion	1. Is the investigation an	exa	am	ple	e of	a	fai	r te	st?	? E	хp	la	in l	οу	ref	err	ing	j to	th	e v	ari	ables
Data	in the experiment.																					
What does your data tell you?																						
data ten you:																						
	2. State which antibacter	ial	ag	ent	Wa	as ·	the	e be	est	an	d١	٧hi	ich	Wá	as ·	the	e w	ors	st ir	ո th	nis	
	experiment.		Ŭ																			
	слренниене.																					

	results support or refute your hypothesis. Justify your response us ata from the results.
	de which antibacterial agent performed the best. Explain your answ
using your	r results to support your answer.

	5. Describe any problems you and the class had with this experiment. Explain what
	changes would you make if you were to repeat it.
	grant gr
Conclusion	
4 sentences	
summarising the experiment.	
1-aim	
2-method 3-results	
4-hypothesis	
	
	