



Chapter [5.37]

Small Vacuum Oven for Polymer Curing**(vacoven2)****(380)****1.0 Equipment Purpose**

- 1.1 The Precision Scientific Vacuum Oven (vacoven2), located adjacent to the SEM coater, is available for baking and curing samples under vacuum.

2.0 Material Controls & Compatibility

- 2.1 vacoven2 can be used for curing and baking samples under vacuum in a range of temperatures from room to 170°C. Members have used it for Polydimethylsiloxane (PDMS), Santovac wafer bonding, curing, offgassing and photoresist bakeout.
- 2.2 Do not place beakers with flammable solvents inside the oven

3.0 Training Procedure & Applicable Documents**3.1 Member Qualify Tool**

- 3.1.1 Timeline (estimated time to completion: 1 day)
- 3.1.2 Train and qualify with any qualified member. All qualified members are superusers on these tools.

4.0 Definitions & Process Terminology

- 4.1 **Curing** - processing meant to toughen or harden polymer material
- 4.2 **Degassing** - removal of unwanted gas from solvents
- 4.3 **Offgassing** - also known as *outgassing*, releasing of trapped gas from material

5.0 Safety

- 5.1 **Burn hazard.** Surfaces inside the oven may burn you. Avoid touching any inside oven surface when hot and wear appropriate hand protection.
- 5.2 Do not place beakers with flammable solvents inside the oven
- 5.3 Use heat resistant red rubber **hot hands** to remove your sample holder
- 5.4 **Researchers must consider the temperature limits of their sample holders.** Some plastics have limited use temperatures; never use a plastic container in a heated oven. Check with staff if you are not sure your sample holder can withstand operating temperature.
- 5.5 The oven may be left unattended when performing well understood, stable processes, during regular MNL Staff hours (8AM – 5PM) if labeled with a blue process identification tag. The oven should be periodically monitored. In the event a process requires oven use outside of regular Staff hours, members are responsible for conducting routine inspections that maintain a safe, well controlled process.

6.0 Process Data

- 6.1 This tool is not monitored for process parameters
- 6.2 Material use history data is collected upon enable and can be found in comment history of the tool in Mercury Client

7.0 Available Processes, Gases, Process Notes

- 7.1 Main processes include heating the sample at room temperature, heating the sample under vacuum and pumping vacuum on the sample at room temperature.
- 7.2 The oven is equipped with (See Figure 11.1)
 - 7.2.1 Temperature control knob
 - 7.2.2 Temperature gauge
 - 7.2.3 Valve to the vacuum pump
 - 7.2.4 Vacuum gauge
 - 7.2.5 Backfill valve

8.0 Equipment Operation

- 8.1 Turn on the Power Switch
- 8.2 Place your sample inside using your own heat resistant cassette or sample holder and take care that your sample does not touch the walls. Close the door.
- 8.3 If vacuum is desired, make sure the backfill valve is closed, then open the pump valve by flipping the black switch up towards the ceiling.
- 8.4 To start the heating process set the temperature control knob to a desired value. Members report a temperature of about 130°C at the setting of 5 to 6. Calibrate the temperature setting for your needs.
- 8.5 After finishing your process, turn the temperature control knob to zero. The oven backfills with N₂. Close the pump valve by flipping the black switch into the horizontal position. Open the backfill valve slowly and let the chamber reach atmosphere (~760Torr). Take precaution not to touch hot inside surfaces when removing your sample. The inside of the door might be hot as well.
- 8.6 Turn off the N₂ flow by closing the backfill valve. Please pump the chamber when you leave the tool.

9.0 Troubleshooting Guidelines


- 9.1 If the tool won't pump down check the door gasket for particles and wipe down with IPA and a Technicloth. Report a fault if condition persists.

10.0 Study Guide


- 10.1 There is no online exam for this tool

11.0 Appendices

11.1 Figures



Front of the oven



Valves at the back