

## Salt Hydrolysis

We expect solutions of substances such as HCl and HNO<sub>2</sub> to be acidic and solutions of NaOH and NH<sub>3</sub> to be basic. However, we may be somewhat surprised at first to discover that aqueous solutions of some salts such as sodium nitrite, NaNO<sub>2</sub>, and potassium acetate, KC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>, are basic, whereas others such as NH<sub>4</sub>Cl and FeCl<sub>3</sub> are acidic. Recall that salts are the products formed in neutralization reactions of acids and bases. For example, when NaOH and HNO<sub>2</sub> (nitrous acid) react, the salt NaNO<sub>2</sub> is formed.

Nearly all salts are strong electrolytes and exist as ions in aqueous solutions. Many ions react with water to produce acidic or basic solutions. The reactions of ions with water are frequently called hydrolysis reactions. We will see that anions such as CN<sup>-</sup> and C<sub>2</sub>H<sub>3</sub>O<sub>2</sub><sup>-</sup> that are the conjugate bases of the weak acids HCN and HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>, respectively, react with water to form OH<sup>-</sup> ions. Cations such as NH<sub>4</sub><sup>+</sup> and Fe<sub>3</sub><sup>+</sup> come from weak bases and react with water to form H<sup>+</sup> ions.

In this lab, you will predict whether various salt solutions will be acidic or basic. You will then test those salts to determine whether your predictions were correct or not.

### Chemicals

0.10 M solutions of the following salts:

AlCl <sub>3</sub>	Na <sub>3</sub> PO <sub>4</sub>	ZnSO <sub>4</sub>
K <sub>3</sub> PO <sub>4</sub>	Ba(NO <sub>3</sub> ) <sub>2</sub>	Na <sub>2</sub> CO <sub>3</sub>
KNO <sub>3</sub>	NaHCO <sub>3</sub>	CuCl <sub>2</sub>
Ca(NO <sub>3</sub> ) <sub>2</sub>	NaCl	NaC <sub>2</sub> H <sub>3</sub> O <sub>2</sub>
Na <sub>2</sub> SO <sub>4</sub>	NH <sub>4</sub> Cl	FeCl <sub>3</sub>

### Equipment

Dropper bottles, Universal pH paper with color codes, scissors

### Procedure

#### **1. Make Predictions**

On your data sheet, predict whether each solution will be acidic, basic, or neutral.

#### **2. Test solutions**

Obtain 4 pieces of universal indicator (orange) pH paper. Cut each pH paper into 4 smaller pieces. Apply a drop of a solution to one of the small pieces of pH paper. Record the pH value indicated by the paper on your data sheet. Repeat for all solutions.

### Disposal

As you test each salt solution, place all used test strips on a piece of paper towel. At the end of the lab, throw away the paper towel and all test strips.

Name: \_\_\_\_\_

### Salt Hydrolysis Lab

#### Predictions:

$\text{AlCl}_3 =$  \_\_\_\_\_

$\text{Na}_3\text{PO}_4 =$  \_\_\_\_\_

$\text{ZnSO}_4 =$  \_\_\_\_\_

$\text{K}_3\text{PO}_4 =$  \_\_\_\_\_

$\text{Ba}(\text{NO}_3)_2 =$  \_\_\_\_\_

$\text{Na}_2\text{CO}_3 =$  \_\_\_\_\_

$\text{KNO}_3 =$  \_\_\_\_\_

$\text{NaHCO}_3 =$  \_\_\_\_\_

$\text{CuCl}_2 =$  \_\_\_\_\_

$\text{Ca}(\text{NO}_3)_2 =$  \_\_\_\_\_

$\text{NaCl} =$  \_\_\_\_\_

$\text{NaC}_2\text{H}_3\text{O}_2 =$  \_\_\_\_\_

$\text{Na}_2\text{SO}_4 =$  \_\_\_\_\_

$\text{NH}_4\text{Cl} =$  \_\_\_\_\_

$\text{FeCl}_3 =$  \_\_\_\_\_

#### Experimental Result:

$\text{AlCl}_3 =$  \_\_\_\_\_

$\text{Na}_3\text{PO}_4 =$  \_\_\_\_\_

$\text{ZnSO}_4 =$  \_\_\_\_\_

$\text{K}_3\text{PO}_4 =$  \_\_\_\_\_

$\text{Ba}(\text{NO}_3)_2 =$  \_\_\_\_\_

$\text{Na}_2\text{CO}_3 =$  \_\_\_\_\_

$\text{KNO}_3 =$  \_\_\_\_\_

$\text{NaHCO}_3 =$  \_\_\_\_\_

$\text{CuCl}_2 =$  \_\_\_\_\_

$\text{Ca}(\text{NO}_3)_2 =$  \_\_\_\_\_

$\text{NaCl} =$  \_\_\_\_\_

$\text{NaC}_2\text{H}_3\text{O}_2 =$  \_\_\_\_\_

$\text{Na}_2\text{SO}_4 =$  \_\_\_\_\_

$\text{NH}_4\text{Cl} =$  \_\_\_\_\_

$\text{FeCl}_3 =$  \_\_\_\_\_

#### Questions:

1) Select one of the acidic salts above and write the hydrolysis reaction for this salt.

2) Select one of the basic salts above and write the hydrolysis reaction for this salt.

3) Calculate the pH of the  $\text{NH}_4\text{Cl}$  solution you used in the lab. How does your answer compare with the experimental result?