


Investigating the Effect of Temperature on the Rate of Hydrolysis of Aspirin

Aim

To investigate the effect of varying the temperature on the rate of hydrolysis of Aspirin using colorimetry.

Background

Aspirin is the trademark for an ester, acetylsalicylic acid, which is a prescription-free painkiller. The active ingredient in Aspirin also prevents infections and damage to blood vessels. The efficiency of any drug depends on its chemical stability. Hydrolysis of the drug can be a major reason for the instability of drug solutions. When Aspirin undergoes hydrolysis, the products are salicylic acid and ethanoic acid.



Salicylic acid forms violet-blue complexes with Fe³⁺ ions. The intensity of the colour depends on the salicylic acid concentration in the sample.

Curriculum Link

Topic IX Rate of Reaction

Topic XV Analytical Chemistry

Chemicals, Apparatus and Equipment

- Aspirin tablet
- Salicylic acid
- Iron(III) nitrate-9-water
- Concentrated nitric acid
- Ethanol (95%)
- Auto pipettes (1 cm³, 5 cm³)
- 25 cm³ beaker x 8
- 100 cm³ beaker x 2
- 100.00 cm³ volumetric flask x 3
- 250.00 cm³ volumetric flask x 1
- 100 cm³ reagent bottle
- Mortar and pestle
- Electronic balance
- Colorimeter

Hazards:

Salicylic acid	
Concentrated nitric acid	 
Ethanol	
Iron(III) nitrate-9-water	 

Safety precautions:

- Wear safety glasses and chemical-resistant gloves
- Carry out the addition of concentrated nitric acid in the fume cupboard

Procedure

Part I: Preparation of a calibration curve

1. Dissolve 0.100g of salicylic acid in a mixture of 10 cm³ ethanol and 10 cm³ of deionised water. Transfer the solution to a 100.00 cm³ volumetric flask. Add deionised water so that the final volume of solution is 100.00 cm³.
2. Transfer 10 cm³ of the solution prepared in Step 1 and dilute it to 100.00 cm³ with 1:1 mixture of 95% ethanol and water in a 100.00 cm³ volumetric flask. The resultant solution is Solution A.
3. Dissolve 4.00 g of Fe(NO₃)₃•9H₂O in 50 cm³ deionised water. Transfer the solution to a 100.00 cm³ volumetric flask. Add 4 cm³ concentrated nitric acid. Add deionised water so that the final volume of solution is 100.00 cm³. The resultant solution is Solution B.
4. Prepare solutions in 25 cm³ beakers as follows:

	1	2	3	4	5	6	7	8
Solution A (cm ³)	0	0.50	1.00	1.50	2.00	2.50	3.00	3.50
Solution B (cm ³)	0.50	0.50	0.50	0.50	0.50	0.50	0.50	0.50
Deionised water (cm ³)	4.00	3.50	3.00	2.50	2.00	1.50	1.00	0.50
Concentration of salicylic acid (mg/dm ³)	0	12.50	25.00	37.50	50.00	62.50	75.00	87.50
Absorbance								

5. Measure the absorbance of solutions using colorimeter (green light of wavelength = 533 nm).
6. Prepare a calibration curve.

Part II: Hydrolysis of Aspirin tablet

1. Weigh one tablet of Aspirin. Grind it into powder using mortar and pestle, and dissolve it in 10 cm³ of 1:1 mixture of ethanol and water.
2. Filter the solution into a 250.00 cm³ volumetric flask, rinse the filter paper with 10 cm³ of 1:1 mixture of 95 % ethanol and water. Add deionized water so that the final volume of solution is 250.00 cm³.
3. Prepare a blank solution of 4.00 cm³ of deionized water and 0.50 cm³ of Solution B.
4. Prepare a solution of 4.00 cm³ of Aspirin solution and 0.50 cm³ of Solution B.
5. Record the absorbance of the Aspirin sample.
6. Repeat the measurement of the absorbance of the Aspirin solution after 15 minutes, 30 minutes, and 45 minutes.

Reference:

The experiment “Comparison of Aspirin Hydrolysis at Ambient and Elevated Temperature” under the project “Hands-on approach to analytical chemistry for vocational schools II: AnalChemVoc II” led by Natasa Gros of University of Ljubljana.