



| | | |
|--|---|---|
|  UNIVERSITY <i>of</i> SAN CARLOS <small>SCIENTIA • VIRTUS • DEVOTIO</small> | Department of Chemical Engineering |  |
|--|---|---|

Experiment Plan
(Form CHEL-1)

| | | |
|---------------------------------------|---|---|
| Laboratory Course | : | ChE2204L |
| Experiment Title | : | Melting Point Determination Using the Thomas-Hoover Melting Point Apparatus |
| Group Code | : | ChE2204LW01 |
| Students' Name & Signature | : | Cabo, Michelle |
| | | Sumarago, Erwin |
| | | Talandron, Rhoel |
| Scheduled Date | : | April 15, 2020 |
| Submission Number | : | 1 |
| Teacher | : | Engr. May V. Tampus |
| Term & Academic Year | : | 2 nd Sem, 2019-2020 |

Teacher's Approval

This is to attest that the students have passed the pre-lab interview and are deemed prepared to conduct the experiment.

**Assessed and Evaluated
By:**

Engr. May V. Tampus

(Signature over printed name)

Date and Time


Objectives of the Experiment

1. To determine the melting point of several compounds using the Thomas-Hoover melting point apparatus; and
2. To evaluate the purity of these compounds.

Information Flow Diagram

Objective 1: To determine the melting point of several compounds using the Thomas-Hoover melting point apparatus.

Calibration of Thomas-Hoover Melting Point Apparatus



Objective 2: To evaluate the purity of these compounds.

$$T_c (C) - T_s (C)$$

Legends:

T_{max} = maximum temperature

T_s = temperature at start of melting

T_c = temperature at complete melting

Work Plan

| Time | Task | Person Responsible |
|---------------|--|--------------------|
| 8:30 – 9:00 | Log-in and Questioning | Cabo |
| | | Sumarago |
| | | Talandron |
| 9:00 – 9:10 | Borrowing, cleaning and preparing of materials to be used in the experiment | Cabo |
| | Preparing the sample. | Talandron |
| | Checking the oil bath and preparation of the Thomas-Hoover melting point apparatus | Sumarago |
| 9:10 – 9:15 | Preparing the new capillary tubes and making sure of proper sealing and are clean and dry. | Talandron |
| | Noting the model number of the melting apparatus. | Cabo |
| | Turning on the stirrer and adjusting until optimum setting is determined. | Sumarago |
| 9:15 – 10: 00 | Turning on the transformer knob setting and starts the timer. | Cabo |
| | Taking the temperature readings. | Sumarago |
| | Recording the temperature readings. | Talandron |
| 10:00 – 10:05 | Plotting the temperature against time for calibration curve. | Sumarago |
| 10:05 – 10:20 | Loading the sample to the capillary tube. | Cabo |
| | Placing the sample into the sample compartments drilled into the cap of the bath. | Sumarago |
| | Determining the heater transformer knob setting and records it. | Talandron |
| 10:20 – 10:21 | Starts the timer. | Cabo |
| 10:21 – 10:35 | Observing the sample until it completely melts. | Sumarago |
| | Records the temperature reading. | Talandron |
| | Turning transformer knob to zero and allow bath to cool. | Sumarago |
| 10:35 – 11:20 | Create a duplicate run. | Cabo |
| | | Sumarago |
| | | Talandron |
| 11:20 – 11:30 | Housekeeping | Cabo |
| | | Sumarago |
| | | Talandron |

Table 1. Thomas-Hoover Melting Point Apparatus Data

| | |
|--------------|--|
| Model Number | |
|--------------|--|

Table 2. Data for the Calibration of Thomas-Hoover Melting Point Apparatus

[illegible]**Table 3.** Transformer Knob Setting

| Knob Setting | Temperature Maximum (°C) |
|--------------|--------------------------|
| 1 | |
| 2 | |
| 3 | |

Table 4. Data for Determination of Melting Point using Thomas-Hoover Melting Point Apparatus

[illegible]

Table 5. Melting Point of Samples

| Sample | Trial | Temperature (°C) | |
|--------|-------|------------------|-------------------|
| | | Started to Melt | Completely Melted |
| 1 | 1 | | |
| | 2 | | |
| 2 | 1 | | |
| | 2 | | |

Observations:



