

General Information for cell culture

A. Preparation

1) Media

- a. Add media (1 pack powder) in 900 ml of dd H₂O
- b. Dissolve NaHCO₃, antibiotic, and NAA with stirring
- c. Adjust volume to 1000 ml dd H₂O
- d. Adjust to pH 7.2 with 1 N HCl
- e. Sterilize by filtration through a 0.2 µm pore size filter under the clean bench
- f. Aliquot into 500 ml sterilized bottle and store 4°C
- g. Add serum and warm in 37°C water bath before use

2) Inactivation of serum

- a. Thaw serum at 37°C in water bath
 - b. Heat at 60°C in water bath for 30 min
 - c. Cool
 - d. Aliquot into 100 ml sterilized bottle and store at -20°C
- serum may also contain substances inhibiting cell proliferation. Some of these maybe artifacts of preparation, e.g., bacterial toxins from contamination prior to filtration; the γ-globulin fraction may contain antibodies cross-reacting with the culture. Heat inactivation removes complement from the serum and reduces the cytotoxic action of immunoglobulin without damaging polypeptide growth factors.

3) Vinblastine solution

- a. Dissolve vinblastine (10 mg/ml) in 100% DMSO and freeze at -20 °C for next use
- b. Add 100 µl vinblastine solution in 1000 ml media (1 mg vinblastine/1000 ml media)

4) Hormone-free serum

- a. Material : activated charcoal, Dextran 70, PBS
- b. Methods
 - i. Prepare the charcoal-dextran solution (250 mg charcoal and 25 mg dextran/50 ml PBS)
 - ii. Add 50 ml charcoal-dextran solution to 1000 ml serum (FBS)
 - iii. Incubate at 55 °C in water bath for 30 min with stirring
 - iv. Allow to cool
 - v. Centrifuge or 10 min at 2000 rpm
 - vi. Collect supernatant
 - vii. Filter sterilize (0.2 µm) (1 filter for 100 ml serum)
 - viii. Aliquot into 100 ml sterilized and store -20 °C

5) Hormone solution

- a. Estradiol for ZR-75-1
 - i. Dissolve estradiol in ethanol and freeze at -20 °C for next use (314.4 µg estradiol/ml ethanol)
 - ii. Add 10 µl of ethanol solution in 100 ml distilled water (3.144 µg estradiol/100 ml distilled water)
 - iii. Filter (0.2 µm) sterilize and store at 4 °C
 - iv. Add 0.5 ml of water solution in 500 ml of RPMI media (31.33 ng estradiol/1000 ml media, 0.1 nM)
- b. Testosterone for LNCaP
 - i. Dissolve testosterone acetate in ethanol and freeze at -20 °C for next use (330.5 µg testosterone acetate/ml ethanol)
 - ii. Add 10 µl of ethanol solution in 100 ml distilled water (3.305 µg testosterone acetate/100 ml distilled water)
 - iii. Filter (0.2 µm) sterilize and store at 4 °C (33.05 ng testosterone acetate/1000 ml media, 0.1 nM)

6) 0.4% Sulforhodamine B (SRB)

- a. Dissolve SRB (0.265 g) in 500 ml acetic acid (1%)

b. Dye content approx, 75%

7) Phosphate-buffered saline (PBS)

- a. Dissolve NaCl (8 gm), KCl (0.2 gm), Na_2HPO_4 (1.15 gm) and KH_2PO_4 (0.2 gm) in 900 ml distilled water with stirring
- b. Adjust volume to 1000 ml with dd H_2O
- c. pH 7.4
- d. autoclave for 30 min on liquid setting
- e. cool liquid
- f. store at room temperature until needed

B. Procedure for freezing cells

- 1) collect cells
- 2) count cell number
- 3) concentrate the cell number (1×10^6 cells/ml media with 10% DMSO)
- 4) transfer 1 ml of cell solution into cryogenic vial
- 5) put vials into Nalgene cryo freezing container which containing isopropyl alcohol
- 6) keep the container at -70°C overnight
- 7) transfer the vials into liquid nitrogen tank

C. Procedure for recovering frozen cells

- 1) add 15 ml of media into a 75 cm^2 tissue culture flask
- 2) thaw the frozen vial in water bath at 37°C
- 3) spray 70% alcohol outside of vial for sterilization
- 4) transfer cells into tissue culture flask
- 5) change the media after 1 day incubation