

CRISPR: sgRNA design and synthesis

Helen Willsey Lab

Updated: August 23, 2023 by Jean Dea

Adapted from: NEB Engen sgRNA synthesis protocol



This protocol was developed with lots of support from the Khokha lab and the NXR. Please cite [this](#), [this](#) and [this](#).

Reagents:

- EnGen sgRNA Synthesis Kit (NEB Product # E3322)
- Zymo RNA Clean & Concentrator kit (Genesee Product # 11-352B)
- Cas9-NLS protein (MacroLabs, UC Berkeley)
- Comes at 6.4mg/ml in 20 mM HEPES-KOH pH 7.5, 150mM KCl, 10% glycerol, 1mM DTT in 10µl aliquots at -80°C. To each aliquot, add 10µl pure water, then aliquot as 1.05µl single-use tubes (1.05µl = 2.2µg)
- Dextran-Alexa555 (Thermo Fisher Product # D34679)

1. DNA oligo design

5' CTAGCTAATACGACTCACTATAG(N)₂₀GTTTtagagCTAGAAATAG 3'
T7 promoter, target, crRNA complement

At least one G must be present at the 5' end of the guide (blue)
sgRNA does NOT contain the PAM sequence.

*** go to crisprscan.org, then click genome tracks at the top, then click *Xenopus tropicalis*. Search your gene. Pick a bright green crispr (bright green = no/few off-targets) that targets an early exon with a high score (>60). Click on it and it will tell you the target sequence and the oligo to order.

2. sgRNA Synthesis (EnGen kit)

☐ At room temperature in the order listed (avoid master mixes):

- | | |
|--|-------------|
| <input type="checkbox"/> RNase-free H ₂ O | 2 µl |
| <input type="checkbox"/> EnGen 2X sgRNA rxn mix | 10 µl |
| <input type="checkbox"/> DNA Oligo (1µM) | 5 µl |
| <input type="checkbox"/> DTT (0.1M) | 1 µl |
| <input type="checkbox"/> Enzyme Mix | <u>2 µl</u> |
- 20 µl reaction

☐ Incubate at 37°C for at least 4 hours.

☐ Add 30µl RNase-free H₂O and 2µl DNase I, mix and incubate at 37°C for 15'.

3. RNA purification (Zymo RNA Clean and Concentrator)

- ☐ Add 2 volumes (100 µl) RNA Binding Buffer to each sample and mix.
- ☐ Add an equal volume (150 µl) of ethanol (95-100%) and mix.
- ☐ Transfer to Zymo Column in Collection Tube, spin 30 seconds. Discard flow-through.
- ☐ Add 400 µl RNA Prep Buffer and centrifuge for 30 seconds. Discard flow-through.
- ☐ Add 700 µl RNA Wash Buffer and centrifuge for 30 seconds. Discard flow-through.
- ☐ Add 400 µl RNA Wash Buffer and centrifuge for 2 minutes.
- ☐ Transfer the column into an RNasefree tube. Add 15 µl RNase-Free H₂O directly to the column matrix and centrifuge for 30 seconds.

- ☐ Nanodrop, dilute to 830 ng/μl in RNase-free H₂O, run 1 μl on gel, and aliquot in single-use volumes (1.5 μl) at -80°C.

4. Injection

Mix:

<input type="checkbox"/> Cas9-NLS protein (2.2 μg)	1.05 μl
<input type="checkbox"/> Dextran-555 (2.5 μg**)	0.5 μl
<input type="checkbox"/> sgRNA (1.2 μg total; 828ng/ul)	<u>1.45 μl</u>
	3 μl

**adjust KCl of dextran so that final injection mix is ~300 mM salt to prevent needle clogging

We mix 20 μl dextran-555 with 20 μl 3.6 M KCl

(3.6 M KCl is in a 15ml Falcon on Helen's bench)

X. tropicalis: Inject 2nl (8 hashmarks) = 1.5 ng Cas9 + 800 pg sgRNA

X. laevis: Inject 4nl (10 hashmarks) = 3ng Cas9 + 1.6ng sgRNA

Target animal pole, as that is where the nucleus is.

CRISPR



PROTOCOL A (Younger than [Stage](#) 44, or if more penetration is needed)

Protocol:

- ☐ 1.
- ☐ 2.
- ☐ 3.

PROTOCOL B ([Stage 44+](#))

Protocol:

- 1) (but make the blocking step overnight, dilute antibodies in 10% CAS in PBT, and then primary is 2hr RT)

CRISPR



PROTOCOL C (Quick cilia [staining](#))

Protocol: