

Titration Lab- Use titration skills to determine what the molarity and the pH are for the unknowns. Show your work. $M_a V_a = M_b V_b$

1. Make 100ml of .1M NaOH using your 1M NaOH. Cover your concentrated 1M NaOH and set aside.
2. Use a 100ml graduated cylinder to place 15ml of unknown acid in an erlenmeyer flask.
3. Rinse out the graduated cylinder very well and place all of your NaOH into it. Use that and/or the 10ml graduated cylinder to keep track of how much NaOH you are adding. Always make starting measurements and ending measurements from here.

Unknown A- ____ ml of the unknown HCl

molarity of A-

pH of A-

Unknown B- ____ ml of unknown HCl

molarity of B-

pH of B-

Unknown C- ____ ml of unknown HCl

molarity of C-

pH of C-

Unknown D- ____ ml of unknown HCl

molarity of D-

pH of D-

Unknown E- ____ ml of unknown HCl

molarity of E-

pH of E-

Unknown F- ____ ml of unknown HCl

molarity of F-

pH of F-