

# XA\_CHO Manual

## Introduction

XA\_CHO is a Graphical User Interface (GUI) designed to streamline the process of generating contrast-detail curves for medical imaging systems. This program assists users throughout the process, from loading images to obtaining final contrast-detail curves.

## Getting Started

To begin using XA\_CHO, follow these steps:

- Launch the XA\_CHO program.
- Load a DICOM multi-frame sequence of the phantom acquired with the desired protocol.

## Features

### GUI Overview


The GUI of XA\_CHO consists of three main phases:

Phase 1: Load Images or Workspace. This phase allows users to load DICOM multi-frame sequences or workspaces.

Phase 2: CHO Contrast-Detail Curve Calculation. Users can perform contrast-detail curve calculations in this phase.

Phase 3: 2-AFC Human Test. This phase facilitates 2-AFC (Two-Alternative Forced Choice) human testing.

### Registration Process



**Figure 1** Leeds TO10 phantom with targets labeled as used in this program.

After loading the DICOM sequence, users can register the phantom using a 6-point affine geometrical transformation. This transformation is calculated based on the indicated six points of the phantom: K1, L1, G1, H1, J1, D1, and D4 targets. These points represent different diameters and contrasts of the phantom targets.


### Manual Indication of Points

Users manually indicate the points on the averaged sequence.

The data from these points are used to calculate the geometrical transformation using a previously calculated parametrical description of the TO10 phantom.

## AFC Human Test


After the registration process, users can perform a 2-AFC (Two-Alternative Forced Choice) human test. In this test, human observers are presented with pairs of images and must indicate whether the signal target is in the left or right image.



**Figure 2** XA\_CHO Analyser GUI.

### Procedure:

1. Presentation of Image Pairs: Pairs of images are presented to human observers.
2. Observer Response: Observers indicate whether the signal target is in the left or right image.
3. Data Collection: The responses from human observers are collected for further analysis.



## Tuning Alpha and Computing CD Curves

Following the 2-AFC test, the alpha value is tuned and set based on the collected observer data. Once the alpha value is set, users can compute a Contrast-Detail (CD) curve using the CHO (Channelized Hotelling Observer) model AUCs (Area Under the Curve).

### Procedure:

1. Alpha Tuning: Based on the collected human observer data, the alpha value is adjusted to optimize the performance of the CHO model.
2. Setting Alpha: Once the optimal alpha value is determined, it is set for further calculations.
3. CD Curve Computation: Using the set alpha value, users can compute the CD curve utilizing the AUCs generated by the CHO model.
4. Storage: The computed CD curves can be stored for future reference or analysis.

## Conclusion

By incorporating the 2-AFC human test and utilizing the CHO model AUCs for CD curve computation, XA\_CHO provides users with a comprehensive tool for evaluating the performance of medical imaging systems. This integrated approach enhances accuracy and reliability in assessing contrast-detail relationships.