

## TZ May 2026 Fischer Lab

### Extract DNA:

- 1) Suspend colonies in 150 mL of distilled water
- 2) Boil in water bath for 10 mins
- 3) Ice bath for 10 mins immediately after boiling
- 4) Centrifuge at 12700 rpm for 2 mins at room temperature
- 5) Transfer 50 mL of the supernatant to an autoclaved tube, store at  $-20^{\circ}\text{C}$

Or use DNA extraction protocol using phenol-chloroform extraction.

### Primers:

- 1) To lyophilized primers from IDT, add listed amount of water to achieve 100  $\mu\text{M}$  conc.
- 2) Before diluting primers in water to achieve a 1  $\mu\text{M}$  primer conc., calculate the total volume required for the number of reactions (+1 extra) + controls.
  - a. For 20 isolates (+1) & 4 controls; 25 total reactions = 250  $\mu\text{L}$  Primer Mix
  - b. For 250  $\mu\text{L}$  Primer Mix @ 1  $\mu\text{M}$  primer concentration =  
2.5  $\mu\text{L}$  of FW primer + 2.5  $\mu\text{L}$  of RV primer + 245  $\mu\text{L}$  ddH<sub>2</sub>O
  - c. For 100  $\mu\text{L}$  Primer Mix @ 1  $\mu\text{M}$  primer concentration =  
1  $\mu\text{L}$  of FW primer + 1  $\mu\text{L}$  of RV primer + 98  $\mu\text{L}$  ddH<sub>2</sub>O
  - d. For multiplex PCR, decrease the amount of water added while keeping the volume of each primer the same.
    - i. Example: 250  $\mu\text{L}$  Primer Mix using **Two** Primer Sets =  
2.5  $\mu\text{L}$  of FW Primer **X** + 2.5  $\mu\text{L}$  of RV Primer **X** +  
2.5  $\mu\text{L}$  of FW Primer **Y** + 2.5  $\mu\text{L}$  of RV Primer **Y** +  
**240**  $\mu\text{L}$  ddH<sub>2</sub>O
- 3) Make 1X Taq MasterMix by diluting your 1  $\mu\text{M}$  Primers 1:1 with 2x Taq MasterMix
  - a. Ex: For 250  $\mu\text{L}$  of 1 $\mu\text{M}$  Primer Mix, add 250  $\mu\text{L}$  of 2x Taq MasterMix
- 4) Transfer 20  $\mu\text{L}$  of 1x MM to each PCR tube
- 5) Add 1.5  $\mu\text{L}$  of Isolate DNA extraction to its labeled PCR tube
- 6) Transfer tubes to thermocycler & follow next steps

### Thermocycler Conditions for PCR using NEB Taq 2x MasterMix:

| Step                                 | Temperature | Time                     |
|--------------------------------------|-------------|--------------------------|
| Initial Denaturation                 | 95 °C       | 30s                      |
| 30 Cycles<br><br>$T_a$ of Primers -> | 95 °C       | 20s                      |
|                                      | 45-68 °C    | 15-60s                   |
|                                      | 68 °C       | 1min per KB of PCR Prod. |
|                                      |             |                          |
| Final Extension                      | 68 °C       | 5mins                    |
| Hold @                               | 4 °C        |                          |

### Thermocycler

- 1) Set thermocycler conditions for primer conditions
- 2) Make sure it runs for 30 cycles
- 3) For Primer Treatment #1:
  - a. Set Cyclic Annealing Temperature ( $T_a$ ) to 50 °C for 30s
  - b. Set Cyclic Extension time to **30s**
- 4) For Primer Treatment #2:
  - a. Set Cyclic Annealing Temperature ( $T_a$ ) to 50 °C for 30s
  - b. Set Cyclic Extension time to **90s**

### PCR Genotyping:

| Gene Amplified                          | Annealing Temp (°C) | Product Lengths    | Identified Strains:                      |
|-----------------------------------------|---------------------|--------------------|------------------------------------------|
| Lanthionine synthetase C family protein | 50                  | 308 bp             | Primer Treatment #1: 1, 5, 8, 30         |
| efflux RND transporter subunit          |                     | 144 bp / 501bp     |                                          |
| mcrC                                    |                     | 211 bp             |                                          |
| aldehyde reductase                      | 50                  | 395 bp             | Primer Treatment #2: 59, 398, 582, Other |
| mtlD                                    |                     | 235 bp             |                                          |
| sucD-lytN                               |                     | ~600 bp / ~1500 bp |                                          |

\*Intergenic space between sucD-lytN adds bases, that's why it is ~ product length

## Gel Electrophoresis – Casting the Gel:

- 1) Make a 1.5% agarose gel:
  - a. To 100 mL of 1x TAE buffer in a flat-bottom glass flask, add 1.5g of agarose.
  - b. Heat in microwave 45s x2, and 30s once. Carefully swirl to mix using a rubber glove between heating times.
- 2) Once melted, let cool for 5 mins, then add 10  $\mu$ L of SYBR Safe Gel Stain and swirl until fully mixed
- 3) Prepare agarose gel casting tray with combs making enough wells for your # of reactions + positive & negative controls + DNA ladders
  - a. Ex: 20 reactions + 4 strain controls + 2 ladders = Minimum 26 wells on the gel
  - b. Make sure your casting tray is level!!**
- 4) Ensure that gel combs are in place and tray is level, carefully pour the mixture into the tray.
- 5) Let gel cool for 15-30 minutes undisturbed before using.

## Loading the Gel:

- 1) Load PCR tube samples with NEB Purple 6X Gel Loading Dye
  - a. 4  $\mu$ L of Gel Loading Dye per 20  $\mu$ L reaction
- 2) Carefully remove gel combs by slowly pulling them straight up and out of the gel
- 3) Remove casting tray and arrange it so the wells face the negative electrode of the gel buffer tank
- 4) Fill gel buffer tank with 1X TAE to fill line, nearing the white electrodes on both ends of the gel (should be over the top of the gel)
- 5) Load first well with 10  $\mu$ L 1kb plus DNA ladder
- 6) Load remaining wells with controls & samples
  - a. Larger wells load 20  $\mu$ L
  - b. Smaller wells load 10 to 15  $\mu$ L
  - c. More sample is better – just don't overfill
- 7) Place gel lid with negative electrode at top of well and positive electrode at bottom of gel on the gel buffer tray; plug in to correct places on power supply
- 8) Run ~100V, checking every 30 minutes for how far the lanes have progressed down the gel
- 9) When lanes have progressed over halfway down the gel, turn off power supply and remove lid
- 10) Carefully transfer gel to Bio-Rad gel imaging machine

11) Adjust contrast/brightness on the machine to maximize band visibility, save to flash-drive, toss gel into biohazard bin and wipe down gel imaging machine platform