

SightStand Smartphone Magnifier Stand:

Last Updated February 18th, 2026.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

SightStand Smartphone Magnifier Stand:

A smartphone stand that magnetically sticks to the back of a smartphone to use the smartphone as a stand magnifier. Available from Maxiaids product code 300302 link. SightStand Information Google Doc link. Scan the QR code below.



**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**



Image description begins. A picture of an iPhone on a SightStand with the Magnifier App open. Size 12 print is shown magnified on the screen. End of image description.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

How to use the SightStand:

- 1. Place the SightStand on the back of the smartphone below the back cameras. Use the provided magnetic sticker if the smartphone or smartphone case does not have the magnetic feature.**
- 2. Pull the two legs out.**
- 3. Open a magnifier app. Examples include the Magnifier App (iPhone only) and the weZoom App (Android only).**
- 4. Use the magnifier app handsfree.**

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

How to use the Clip-On Focus Lens:

For smartphones that don't focus with just the SightStand, users must attach the included clip-on focus lens.

Take both covers off of the lens and twist the lens clockwise into the clip. Open a magnifier app and test each of the back cameras to check which will be used by covering it with a finger or something similar.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

The camera that when covered blocks the view is the one that the clip-on focus lens will go on. Lastly, squeeze open the clip-on focus lens and slide it on to the smartphone either from the top side or the left side.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**




Image descriptions begin. Two pictures showing what the clip-on focus lens for the SightStand looks like with and without the covers on. End of image descriptions.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**



Image descriptions begin. Two pictures showing what the clip-on focus lens for the SightStand looks like put together and on an iPhone. End of image descriptions.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

The Three SightStand Positions:

- **The Stand Above Position.** Open both legs as wide as they go. Lay the legs on a flat surface.
- **The Standing Position.** Open both legs as wide as they go. Lay the bottom side of the SightStand on a flat surface with the legs behind it.
- **The Angled Position.** Open both legs as wide as they go. Then push both legs together until they are touching. Then place the SightStand on its side on a flat surface.




Image description begins. A picture of an iPhone on a SightStand with the Magnifier App open. Size 12 print is shown magnified on the screen. End of image description.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**



Image description begins. An iPhone on a SightStand in the upright position. End of image description.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**





Image description begins. A picture of the SightStand. End of image description.

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**



**Image description begins. A picture of the
SightStand. End of image description.**

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**

Online Resources for Information About the SightStand:

- **SightStand website link.**
- **SightStand Help Videos website page link.**
- **Maxiaids SightStand website page link.**
- **SightStand Amazon page link.**
- **SightStand YouTube channel link.**
- **SightStand review YouTube video from The Blind Life link.**
- **SightStand Updates YouTube video from The Blind Life link.**

**Created by Elijah Leightcap, Vision Rehabilitation Therapist,
Bureau of Blindness and Visual Services, Philadelphia, vrtinfo.com**