







PowerMic Mobile - Day to Day Use

PowerMic Mobile (PMM) turns your iPhone into a secure, wireless microphone for use with Dragon Medical One. Below are the steps on how to open PMM for use and how to close it.



Turn On & Off

Access the Dragon's Login page: - Click the Dragon Bar's black flame icon. A menu will appear. - Click log off - Wait a few seconds and the Login Page will automatically appear. - Change microphone to: Nuance PowerMic Mobile	
The Dragon bar will appear on the Epic screen with a new telephone icon. The icon color will change: red/Orange = Microphone is Off green = Microphone is On	
Locate the PowerMic Mobile icon on your phone and double click on it. This will open up the application.	
Log into the phone's PowerMic Mobile application. The PMM User Name is equal to your Epic and Dragon's Userid values.	

Start recording by one tap to the PMM microphone button. Tap again to stop recording. No need to press and hold the button. Note: - Do <u>not</u> press the DAX CoPilot button. This is not a functional button at this time. - PowerMic Mobile supports the use of wired headsets, Bluetooth headsets and earbuds.	Record screen - microphone off  Record screen - microphone on 
Note: If you see an orange icon on the Dragon tool bar...that means the phone needs to be sync'd. Tap the phone microphone button again. The Dragon bar phone icon will change to the correct colors: red = Dragon cannot hear you green = Dragon can hear you.	 
When you are finished with using PowerMic Mobile on the phone: 'Swipe' the application to close.	
TIPs: You can set a notification sound for when you start & stop recording. Steps on setting this: - Click the PowerMic Mobile menu (3 dots) - Click Settings - Select your sound settings - Tap the back button on your device Click the <u>Help</u> feature to view some detailed information on using PMM.	