



## Instruction Sheet

### #303240 - Vibro-Acoustic Floor Mat

**Take time to familiarise yourself with the use and maintenance of this equipment before using. Please file this instruction sheet in an accessible location**

Materials List	Tools Needed
1 x Floor Mat	N/A
Coaxial Cable	

## Dimensions

150 x 90 x 15cm

## Installation

The Floor mat comes delivered completely assembled. In order to make the mat vibrate the Rocker must be connected to a Southpaw Radio Rover or Stereo. Use with any other stereo equipment will void the warranty.

The Floor mat comes with a coaxial cable with RCA jacks at each end. Locate the speaker connection on the bottom of the mat and plug one end of the coaxial cable wire into the RCA socket. Connect the other end of the wire to the RCA socket on a Radio Rover or stereo. If the stereo you are using does not have an RCA socket you can connect the wire to Red/Black posts or connectors by cutting the RCA jack off the speaker wire. Strip 1/2" of insulation to expose the wire. Connect the wire with the black line to the black post or connector. Connect the wire with no black line to the red post or connector.

## Construction

This item is constructed of piece high density foam & a high tear strength, non toxic, phthalate free vinyl.

## **Maintenance**

Checks should be carried out to observe wear points on a monthly basis and recorded on a checklist.

For cleaning, please check with your infection control department and with the client for any allergies before using these procedures.

To clean, wipe with a damp cloth and a mild surface cleaning solution and either let air dry or wipe with a dry cloth.

Cover can be removed for cleaning.

## **Safety**

The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

Products supplied by Southpaw, where necessary, have been designed, manufactured and checked to ensure that they meet EU safety, health and environmental requirements.

Where legal CE marking is required to adhere to the relevant directives or standards, Southpaw certify that these products are CE marked as a declaration of our responsibility and conformity with all legal requirements to achieve CE marking. We therefore can confirm that these products are valid for resale throughout the EU.