

**Annexure-II**

**Calibration Information Sheet for Thermostatic Column Compartment**

<div>Company Logo Here</div> <div><b>XX PHARMACEUTICALS LIMITED</b></div> <div>117 Adams Street, Brooklyn, NY 11201, USA</div>
--

**Calibration Information Sheet for Thermostatic Column Compartment**

Instrument Name:

Instrument ID No.:

Location	Frequency
Instrumental lab	Once in every six month

**Column Temperature Accuracy**

Flow Rate : 1.0 ml/min				
Set Temperature	Observe Temperature	Accuracy	Limit	Result
80.0°C			≤ 3°C	
40.0°C				

**Column Temperature Stability**

Flow Rate : 1.0 ml/min							
Set Temperature	Observe Temperature	Average	Low	High	Stability	Limit	Result
80.0°C						≤ 1°C	
40.0°C							

Calibration Date :	Next Calibration Date :
--------------------	-------------------------

Done By :	Checked By :
Name :	Name :
Designation :	Designation :
Date :	Date :