

Annexure-III

Calibration Information Sheet for Electric Analytical Balance

Company Logo Here XX PHARMACEUTICALS LIMITED 117 Adams Street, Brooklyn, NY 11201, USA

Calibration Information Sheet for Electric Analytical Balance

Instrument Name:

Instrument ID:

Location: Analytical lab	Frequency: Once in every six months
---------------------------------	--

1. Eccentricity:

Date	Test weight (g)	Position	Measurement (g)	Eccentricity Error (g)	Tolerance(g)	Remarks
	50	Center			±0.0003	
	50	Left rear				
	50	Right Rear				
	50	Right Front				
	50	Left Front				
	50	Center				

2. Linearity:

Date	Test weight (g)	Display Value(g)	Deviation(g)	Linearity Error (g)	Tolerance (g)	Remarks
	0				±0.0002	
	10					
	20					
	50					

3. Sensitivity:

Date	Test weight (g)	Display Value(g)	Difference(g)	Tolerance (g)	Remarks
	100			±0.0008	

4. Repeatability (100g):

Date	Empty (g)	Load (g)	Difference (g)	Standard Deviation (g)	Tolerance(g)	Remarks
	0.00000				±0.0001	
	0.00000					
	0.00000					
	0.00000					
	0.00000					
	0.00000					
	0.00000					
	0.00000					
	0.00000					
	0.00000					

Done By :	Checked By :
Name :	Name :
Designation :	Designation :
Date :	Date :