

Annexure-I

Calibration Information Sheet for Fume Hood

Company Logo Here XX PHARMACEUTICALS LIMITED <small>117 Adams Street, Brooklyn, NY 11201, USA</small>

Calibration Information Sheet for Fume Hood

Equipment Name:

Equipment ID.:

Location			Frequency					
Analytical lab			Once in every six months					
Analytical Data (Air flow)								
Flow Meter Reference No.:								
Sl. No.	Position	Flow meter reading (m/sec)	Average reading (m/sec)	Acceptance criteria (Face velocity)	Calibration Pass/Fail			
1.	Left			Not less than 0.5 m/sec from outside to inside. (When the window sash is open 30 cm from the bottom)				
2.	Middle							
3.	Right							

Done By :	Checked By :
Name :	Name :
Designation :	Designation :
Date :	Date :