

IHE Change Proposal

Tracking information:

IHE Domain	Devices (DEV)
Change Proposal ID:	CP-DEV-xxx
Change Proposal Status:	Draft
Date of last update:	2025-03-24
Person assigned:	Silas Morris

Change Proposal Summary information:

OBX MDC_PUMP_VOL_DELIV_DISPLAYED	
Submitter's Name(s) and e-mail address(es):	Silas Morris silas@epic.com
Submission Date:	2025-03-xx
Integration Profile(s) affected:	Infusion Pump Event Communication (IPEC)
Actor(s) affected:	DOR
IHE Technical Framework or Supplement modified:	IHE_PCDTF Vol 2, dated 2024-11-04
Volume(s) and Section(s) affected:	Changes to sections Appendix M, M.1.2
Rationale for Change: If the clinical user ever needs to compare the infused volume on the DOR with the infused volume reported to the EHR, they may be unable to effectively do so if the DOR is reporting volume in MDC_VOL_FLUID_DELIV that doesn't match what is displayed on the device.	

In *Table M.1.2.1-1 – Mapping of Infusion Pump Event Parameters to MDC Codes*, add the following content on page 231, at line 5518:

New content

Parameter	MDC Code
Infused Volume Displayed	MDCX_PUMP_VOL_DELIV_DISPLAYED