

AAC App Feature Matching Comparison

Feature	LAMP WFL	TD Core (Previously Snap+Core First)	TouchChat	Proloquo2go	Proloquo4text	GoTalk Now
Symbol Set	Minspeak (Unity)	Mayer-Johnson PCS/Boardmaker	Symbol Stix	Symbol Stix	None. Text-based only	Mayer-Johnson PCS/Boardmaker
Grid Size	84 icons/page	Up to 80 icons/page -Grid Sizes vary from 1x1 to 8x10	- <i>Word Power</i> - 20/25/42/48/60/80/ 108/120/140 icons/page Alternative Options: - <i>MultiChat</i> 15- 15 icons/page - <i>My QuickChat</i> - 12 icons/page - <i>4 Basic</i> - 4 icons/page	Up to 144 icons/page	N/A, single screen payout	1, 2, 4, 9, 16, 25, 36 icons/page
Organization	Semantic Compaction	Category-based	Category-based	Category-based	N/A, text-based	Category-based
Access Methods	-Direct selection- Touch -Auditory Scanning/Switch Access	-Direct selection- Touch -Auditory Scanning/Switch Access -Eye gaze (Surface Pro)	-Direct selection- Touch -Auditory Scanning/Switch Access -Head tracking (iPad Pro)	-Direct selection- Touch -Auditory Scanning/Switch Access	-Direct selection- Touch -compatible with accessible keyboards and switches	-Direct Selection- Touch -Auditory Scanning/Switch Access
Robust Core Vocabulary?	Yes	Yes	Yes	Basic, Intermediate, and Advanced levels	Yes	Basic Core Vocabulary

High Contrast Icons?	Yes-manually customize	Yes	Yes-manually customize	Yes-manually customize	No	Yes-manually customize
Feature	LAMP WFL	TD Core (Previously Snap+Core First)	TouchChat	Proloquo2go	Proloquo4text	GoTalk Now
Keyboard function?	Yes	Yes	Yes, with predictive text	Yes	Yes, with predictive text	Yes, with predictive text
Visual Scene Display?	No	No	Yes	No.	No	Yes
Search Option to Assist with Finding Icons	Yes-Word Finder Tool	Yes, magnifying glass at top left of screen	Yes, within the keyboard page	Yes, magnifying glass at bottom of page	No	No
Editing Restriction Settings	Yes, editing menu passcode within app	Yes, editing menu passcode within app	Yes, editing menu passcode within app	Yes, passcode within app and removal of editing icons within iPad Settings App	Yes, editing menu passcode within app	Yes, within Settings on initial page
Additional Touch Activation Settings	-Dwell time settings for activations -Release time settings for restricting repeated hits	-Variety of touch settings with dwell time, connecting multiple taps, delay between selections	-Dwell time settings for activations -Release time settings for restricting repeated hits	-Various Repeated tap customization settings -Repeat Delay option	No	-Debounce Setting: Allows for specified delay between hits
Ability to Increase Button Margins?	Yes, within Settings and Buttons menu	Yes, within Settings>Page Set>Style>Pages	Yes, within the Settings and Buttons menu	Yes, within Settings>Appearance> View	n/a	No
Keyguard Option?	Yes, consult AT Dept	Yes, consult AT Dept	Yes, consult AT Dept	Yes, consult AT Dept	No	Yes, consult AT Dept

Data Tracking Capability?	Yes, consult AT Dept	Yes, with distinction of modeling for user activations	Yes, consult AT Dept	No	No	No
Feature	LAMP WFL	TD Core (Previously Snap+Core First)	TouchChat	Proloquo2go	Proloquo4text	GoTalk Now
Other	No more than 3 button hits to access every icon in the app	Navigation bar allows for quick changes between “core vocab” page, “quick fires” page, fringe vocab/topic pages, and keyboard.	Prioritizes syntax development through linked icons/pages to build sentences and supports literacy development by using a combination of symbols and word-only icons.		Can create and store personalized phrases	This app is fully customized for the user and does not have preloaded pages.