




# Guide for trialing ServiceM8 app updates

---

## How to install the update

- Tap the link provided in the email on the device you want to install the update.
- Install the Apple Testflight app
- Tap the link provided in the email again, and it will open the TestFlight app.
- You'll be presented with the option to accept ServiceM8.



- Once you accept, you'll be able to install the new version of the ServiceM8 app.

## What to do if you have issues with the update

If you're experiencing an issue with the app that is disrupting your work, you can revert back to the store version by deleting the new version and re-installing ServiceM8 from the App Store.

Always make sure your app has fully synced all changes before deleting it, by pulling down the activity feed screen.