

NDH Hospital, NTEC
Contour plus (BGA) Internal QC Record

Clinic / Unit / Ward: OT _ Machine Number of
BGA: 1 / 2 / 3 / 4 Hi / Lo Control Solutions
Expiry Date: Opened on:
To be discarded on:

Perform QC testing at least once daily / before use, repeat under the following conditions:

- A new bottle / pack of strip is opened
- QC testing result is out of the acceptable range.
- Patient's CBG result is questionable.
- BGA is dropped.
- Accuracy of blood ketone result is questioned.
- Battery of BGA has been changed.

Please change a new form when changing a new BGA / control solution

[illegible]

Checking Record for Contour Plus Blood Glucose Meter

Machine No.	EAM	Serial Number	Period of use
1	2568378	7887E212890 4	July- Dec
2	2568380	7887E212890 3	July - Dec
3	2568381	7887E212890 2	Jan - Jun
4	2568379	7887E212890 1	Jan - Jun