



SPECIMEN REQUIREMENTS

For Respiratory Tract Infection PCR Testing

CLIA Certified and CAP Accredited Laboratory

Respiratory Tract Infection (RTI) PCR Testing	
<i>Use</i>	<i>Respiratory Tract Infection PCR Testing</i>
<i>Specimen Type</i>	Anterior nasal swab, Nasopharyngeal swab, oropharyngeal (throat) swab
<i>Container</i>	Universal transport medium (UTM)
<i>Volume</i>	3 ml
<i>Minimal Volume</i>	1 ml
<i>Handling</i>	Keep the samples at 2-8°C after collection. Ship the specimen with a cold pack (2-8°C) to the laboratory on the day of collection for next-day delivery.

Additional Instructions

Please call (877) 888-2973 or email info@virantdx.com for questions and shipping assistance.
Please visit <https://virantdx.com/testing-solutions/> for our complete test menu and detailed test information.

Virant Diagnostics, Inc.

11002 Veirs Mill Rd, #404, Wheaton, MD 20902
Website: virantdx.com Email: info@virantdx.com
Phone: (877) 888-2973 Fax: 888-713-3456

Document ID: VDX-SPE-329-v1
Last modified: 10/06/2023