

Annexure-II

Chemical Calibration Information Sheet for Dissolution Tester

Company Logo Here XX PHARMACEUTICALS LIMITED 117 Adams Street, Brooklyn, NY 11201, USA

Chemical Calibration Information Sheet for Dissolution Tester

Instrument Name:

Instrument ID.:

Location: Analytical lab					Frequency: Once in a year				
Prednisone Std. Wt. (about 20 mg): _____mg					$\% \text{ Dissolved} = \frac{\text{Spl. Abs} \times \text{Std. Wt} \times 5 \times 500 \times \text{Std. Potency} \times 100}{\text{Std. Abs} \times 100 \times 100 \times 10 \times 100}$				
Standard Absorbance at 242 nm: _____									
(Apparatus 1 - Basket)	Tab 1	Tab 2	Tab 3	Tab 4	Tab 5	Tab 6	Tab 7	Tab 8	Limit
Abs. of sample									
% Dissolved									
(Apparatus 2 - Paddle)	Tab 1	Tab 2	Tab 3	Tab 4	Tab 5	Tab 6	Tab 7	Tab 8	Limit
Abs. of sample									
% Dissolved									

Remarks:

Done By	:		Checked By	:	
Name	:		Name	:	
Designation	:		Designation	:	
Date	:		Date	:	