

Lumbar Laminectomy (for removal of tumor)

Indication:	Tumor excision
Position:	Prone with both hands keep in swimming position
Anae:	Remifentanyl x 2 Plain ETT +/- reinforcement tube NS with plain set, SCD, Foley insertion
LA infiltration Please LABEL	<ul style="list-style-type: none"> ✓ Adrenaline (1:10,000) 10ml /amp x 1 ✓ 10 ml Normal Saline x 2 amp ✓ 5mg/ml Chirocaine (0.5%) 10 ml/vial x 1 ✓ 3ml syringe x 1, 20ml syringe x 1, 21G needle x 1 ✓ Skin marker x 1, L sterile stripe x labelling
Instruments	N001, N004, N002b, N106 or N105 N301 N503, N504, N509 N801, N804, N807 +/- N201 for tumor removal
Supple Instruments:	Bipolar forceps (Aesculap) Spinal process bone cutter Tumor holding forceps 1mm, 4mm 1@ Lamp handle x 2 USG probe
Consumables	mini raytec gauze, self adhesive drape 75 x75 4pcs Ioban 6640, Blade No.20, No.11 1@ Bulb syringe x 1, 20ml syringe on blunt needle x 2 set Bone wax, Surgicel S x 1 Pettis M, S 1@, Mic drape, Telfa (for specimen sending frozen section)
Sutures	2/0 Vicryl (muscle closure) Skin stapler (Skin closure)

Operating procedures:

- 1 GA □ 2. Positioning □ 3. Skin preparation
- 2 Draping
 - 2.1 Sterile X-ray cover x 2 to C-arm
 - 2.2 self adhesive drape 75 x75 x 4pcs
 - 2.3 Ioban 6640
 - 2.4 Large disposable drape sheet x 1 to foot end, x 1 to head side
 - 2.5 Small disposable drape sheet to left side and right side. Thus, making a square disinfected area for operation.
- 3 Incision
- 4 Closing