

Company Logo Here

XX PHARMACEUTICALS LIMITED

117 Adams Street, Brooklyn, NY 11201, USA

**Quality Control Department
Microbiology Section
Moist Heat Sterilizer Validation Report**

Section : Date of Test :
Machine ID : Reporting Date :
Sterilizer Condition : Sterilization Temperature :
Exposure Time : Spore Used :
Test Frequency : Next Due Date :
Incubation : 50°C to 55°C/7 days
Temperature/Duration

Sl. No. of spore	Location of spore exposed in the sterilizer	Observation						
		Day-1	Day-2	Day-3	Day-4	Day-5	Day-6	Day-7
1	Top Front inside Load (Non Sterile Side)							
2	Middle Front (Non Sterile Side)							
3	Bottom Front inside Load (Non Sterile Side)							
4	Bottom Back inside Load (Sterile Side)							
5	Middle Back (Sterile Side)							
6	Top Back inside Load (Sterile Side)							
7	Drain Sensor							
8	Chamber Center							
Test Control	Positive Control							
	Negative Control							

“—” sign No growth, “+” sign found growth.

Remarks: Moist Heat Sterilizer is valid/not valid/test to be repeated.

Tested by: _____

Checked by: _____

Approved by : _____

Head of Quality Assurance

--

www.pharmaguideline.co.

Prepared by :	Checked by :
---------------	--------------