




LZ Veto + Offline Data Quality

Ewan Fraser

1

LUX-ZEPLIN!

- Dual Phase LXe TPC.
- Looking for WIMPs. Dark Matter Candidate.
- Surrounded by 2 vetoes. LXe Skin and an Gd-LS Outer Detector.



Offline Data Quality - PREM - Physics Readiness Monitor

Motivation:


- □ Validate that LZ's data Runs are ready for their respective physics analysis.
- □ Verify the stability of the LZ detectors.
- □ Identity and catalogue any defects in the data.


Components:

- □ Automatically run, data analysis modules (ALPACA).
- Data Quality Leads write and coordinate these analysis. Outer Detector Me.
- □ Shifters monitoring of PREM.
- □ Website to display the results of these analysis.
- Database store associated metadata.

ALPACA - LZ's Analysis Package

- Modify ALPACA to output JSON metadata files.
- Analyse time varying parameters.
 - Across multiple Runs.
 - Within a Run.
- Flag significant changes in time varying parameters.






PREMweb

PREM Database:

- Store JSON Metadata files for each run.
- □ Semi-structured, MongoDB

PREM Website:

- □ Frontend: Bootstrap CSS and Javascript.
- Backend: NodeJS App.
- □ Express + Mongoose to interact with DB.



PREM Stack deployed at NERSC

	State 众	Name ᅌ		Image 💠
Namespace: premweb				
	Active	app 🚷 80/http	 Contains frontend and backend. 	efraser/premweb-nodejs 1 Pod / Created a month ago / Pod Restarts: 0
	Active	db 🍪 32457/tcp	←MongoDB.	monqo:4 1 Pod / Created 2 months ago / Pod Restarts: 0
	Active	web 💩 80/http, 443/https	-Hosts site and 2-factor authentication.	IbInest/nginx:1.18.0_1 1 Pod / Created a month ago / Pod Restarts: 0

PREM Stack deployed on SPIN2, viewed in Rancher2.

SPIN2:

- □ NERSC's container-based platform for science gateways.
- Deploy each project using a Docker image.
- Creates a Kubernetes cluster.

Photos! Outer Detector installation



Alice installing the OCS fibres to monitor scintillator and acrylic degradation.



Harvey working on OCS and OD PMT installation.


OD PMT QA: Analysis

Motivation:

- Provide QA for the OD PMTs.
- $\Box \quad Calibrate \longrightarrow Improve our model of the OD.$

Simulations:

- □ SPE: 1k photons, 100k events.
- □ Afterpulsing: 10k photons, 500k events.




OD PMT QA: SPE


- Measure the response of the OD PMTs to a single photoelectron (SPE).
- Spectra in LZ are measured in photons detected (phd).
- Calibrating this response is the first step towards energy calibration.

Analysis Note - Qing He, Daya Bay

$$S_{ideal}(x) = P(n;\mu) \otimes G_n(x) = \sum_{n=0}^{\infty} \frac{\mu^n e^{-\mu}}{n!} \frac{1}{\sigma_1 \sqrt{2n\pi}} \exp(-\frac{(x-nQ_1)^2}{2n\sigma_1^2}).$$



- Photoelectron ionizes a molecule when travelling between photocathode and first dynode.
- The more molecules the more afterpulsing.
- Measure of the quality of the vacuum in a PMT.
- Has it degraded since production in Korea?



Neutron Selection in the LZ Veto

- Begun to look at extending the Neutron Selection analysis for the LZ Veto.
- Improve our understanding of the dominant background in the Veto, gammas from cavern rock.
- Simulating the U238, Th232 and K40 decay chains in the cavern Rock. CPU intensive, not many make it through the water shield.
- Use event biasing to improve the efficiency of the simulation.





Diagram of surfaces in LZ used for event biasing, Andrew N - Sheffield.

Summary

- Developed LZ's Offline Data Quality tool PREM.
- Helped coordinate Data Quality analysis group and the various analysis modules produced.
- Responsible for the Data Quality analysis in the Outer Detector.
- OD PMT QA is about to begin!
- Developed an SPE analysis and Afterpulsing analysis for this effort.
- Begun work on looking at the neutron selection in the LZ Veto.
- Simulating the main veto backgrounds using an Event Biasing technique.

- Analysis: Andrew S Oxford LZ
- Electron Lifetime: Measure of LXe purity Important correction for analysis



Current Development for PREM

Sam Woodford - Liverpool:Anna D - UCL:Automation of analysis.Creating a data structure.Database Queries!Handle calibration values for all PMTs.