

'Flask Mk I' Prosthetic

Based off the 'Lightweight' prosthetic, the 'Flask Mk I' has a cap on the shoulder, which can be removed in order to fill the arm with a liquid of the user's choosing. This liquid is stored in a stain-resistant, flexible container that can hold 16 fl oz on average. This container is refrigerated, causing the exterior of the arm to be cool to the touch. In order to drink from the arm, one need only turn the thumb counter clockwise, allowing it to open up to reveal a nozzle where one may drink the contents from. Notable is that the 'Flask Mk I' can only lift, on average, 80lbs.
