

CT Thorax, Abdomen, and Pelvis w IV contrast **Routine**

1. Thin w Thorax: 1.25 mm standard soft tissue window axial
2. Thick w Thorax: 5 mm standard soft tissue window axial
3. Thin w Thorax: 1.25 mm bone plus lung window axial
4. Thin w Abd/Pel: 1.25 mm standard soft tissue window axial
5. Thick w Abd/Pel: 5 mm standard soft tissue window axial
6. Coronal w Thorax: 3 mm standard soft tissue window
7. Sagittal w Thorax: 3 mm standard soft tissue window
8. Coronal w Abd/Pel: 3 mm standard soft tissue window
9. Sagittal w Abd/Pel: 3 mm standard soft tissue window