

PLEASE CALL DR. JACQUELINE ANGEL (901-481-6787) WITH QUESTIONS. THANKS!

CT Neck With Head Tilt

1. 0.625 mm standard algorithm in soft tissue window (FULL NECK scan) - axial
2. 2.5 mm standard algorithm in soft tissue window (FULL NECK scan) – axial
3. 0.625 mm boneplus algorithm in bone window (FULL NECK scan) – axial
4. 2.5 mm standard algorithm in soft tissue window (FULL NECK scan) – coronal
5. 2.5 mm standard algorithm in soft tissue window(FULL NECK scan) – sagittal
6. 0.625 mm stad algorithm in soft tissue window (HEAD TILT DOWN scan) – axial
7. 0.625 mm standard algorithm in soft tissue window (HEAD TILT UP scan) – axial