
Acute Stroke Brain WO Protocol

Scan notes:


Revised 01/18/2023
Charge as: Brain WO IV CONTRAST
Scanner: 1.5T or 3T
Coil: Head Coil/Flex Coil (flex coil may be needed if anterior head coil cannot clamp)

Plane	Weighting	Slice	Gap	FAT SAT	Scan Range
SAG	T1 MPRAGE 3D	1mm	0.5mm	None	Whole brain/parallel to midline of brain. NO REFORMATS
AXIAL	Diffusion	4mm	0mm	None	Whole brain/parallel to AC/PC line
AXIAL	T2 Flair	4mm	0mm	Strong	Whole brain/parallel to AC/PC line
AXIAL	T2 GRE	4mm	0mm	None	Whole brain/parallel to AC/PC line

SAG Image Planning:



AXIAL Image Planning:



COR Image Planning:

