

*ETN	55	55.00%
Al(3µm/)	12	12.00%
Mg/Al(44µm/325mesh)	18	18.00%
Soy Lecithin(liquid)	0.3	0.30%
HTPB Resin	11.687	11.70%
Isodecyl Pelargonate Plasticizer	0.735	0.74%
Modified MDI Isocyanate Curative	1.838	1.84%
CAO-5 Antioxidant	0.441	0.44%
<p><b>* If you don't use the positive OB% of ETN then you will need to adjust the formula as follows: 45% PETN or RDX and 10% Ammonium Perchlorate.</b></p>		